

Fresenius Medical Care North America Service Bulletin

Equipment: All 2008® Series Hemodialysis Machines
 Bulletin: 22-FRS-002 Rev A
 Subject: Power Supply Heater Cable Wire Connections

1.0 PURPOSE

To inform all Field Service and Technical Support personnel as well as Fresenius Medical Care-trained customers regarding the power supply heater cable wire connections to the distribution board heater block when replacing the cable.


2.0 SCOPE

Danger: Extreme shock hazard. Make sure the power cord has been removed from the wall outlet and the heater switch and main power switch are in their off positions on the rear of the power supply.

Caution: Care must be taken when inserting wires into the Distribution Board Heater Block. Observe the following:

- Examine each wire and remove any stray wires not captured by the ferrule.
- Insert the wire/ferrule into the Distribution Board Heater Block just far enough so no insulation will be crimp when the screw is tightened.
- Before tightening each screw verify again that there are no stray wire strands.
- Securely tighten each screw.

When replacing the power supply heater cable use the illustrations below to assist in correct wire color placement.



If there are any questions regarding this bulletin, contact Fresenius Medical Care Technical Support at 800-227-2572.

© 2022 Fresenius Medical Care. All Rights Reserved.