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CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

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Instructions for Use 650g/900g *bibag*® Bicarbonate Concentrate

2360C2011 SSA

Powder for production of liquid sodium bicarbonate concentrate for bicarbonate hemodialysis

INDICATIONS FOR USE: The *bibag* system is indicated for use in patients undergoing extracorporeal bicarbonate hemodialysis for acute and chronic renal failure. The *bibag* system is intended to be used as one component in the preparation of dialysate according to a physician's prescription in a 3-stream proportioning hemodialysis machine equipped with the *bibag* module.

CONTRAINDICATIONS: There are no absolute contraindications for bicarbonate dialysis. *bibag* bicarbonate concentrate must be used only on the direction of a physician who has evaluated all the pertinent features of the product in relation to the individual patient.

REQUIREMENTS:

- For use with hemodialysis machines with three-stream proportioning bicarbonate systems when set to proportion 1 part acid to 44 parts bicarbonate base and water.
- For use only with systems equipped with the *bibag* module and DIASAFE® *plus* filter.
- For use only with 45X acid concentrates.
- For use only with water that meets current ANSI/AAMI RD62 or ISO 13959 requirements.

WARNING:

- Failure to follow these Instructions for Use may cause patient injury or death.
- Use without any equipment listed in the REQUIREMENTS section or with other equipment than listed in the REQUIREMENTS section may cause patient injury or death.
- If required by the Hemodialysis Machine Instructions for Use, test the conductivity and approximate pH of the dialysate with an independent device before beginning treatment. Use of a dialysate with incorrect conductivity or pH can cause serious injury or death.
- The dissolved *bibag* bicarbonate concentrate must be used within 24 hours of connecting to the dialysis machine. Bacterial growth may occur when using bicarbonate concentrate.
- FOR USE ONLY WITH A HEMODIALYSIS MACHINE EQUIPPED WITH THE BIBAG MODULE. NOT FOR USE BY ANY OTHER ROUTE.

Refer to the Hemodialysis Machine Instructions for Use for a complete list of Warnings associated with the use of *bibag* in hemodialysis treatment.

CAUTION:

- Do not use if the product is open or damaged. 
- Replace a leaking *bibag* immediately. Spills can cause slips and falls; clean up spills immediately.
- Product must be used before the expiration date on the product's label.

ADVERSE REACTIONS: Some undesirable effects related to dialysis can occur, for instance, hypotension, nausea, headache and cramps.

DIRECTIONS FOR USE:

Note: *bibag* bicarbonate concentrate is sufficient for a dialysis session of 6 hours minimum (650g) or 8 hours minimum (900g) at a dialysate flow rate of 500 mL/min and a dialysate bicarbonate setting of 32 mEq/L (mmol/L). Refer to the hemodialysis machine Instructions for Use for treatment times at other dialysate flow rates and bicarbonate settings.

1. Remove the plastic seal from underneath the water and bicarbonate nozzles of the *bibag*. Then open the *bibag* door on the machine by lifting up on the handle.
2. With the white *bibag* handle facing outward, hang the *bibag* on the *bibag* connector nozzles. Close the door, making sure it latches firmly in place.
3. Refer to the detailed directions for use provided in the hemodialysis machine Instructions for Use.
4. At the end of treatment, use the "Empty *bibag*" function of the dialysis machine (or drain to a convenient sink) to empty the bag for easier disposal.
5. Lift up on the *bibag* door handle to open the *bibag* door. Remove the *bibag* and dispose of it and any unused powder according to dialysis unit procedures.
 - Dispose of the bag and unused powder between each treatment.
 - In some cases the bicarbonate powder may contain particles that have been slightly discolored after exposure to water. These discolorations originate from the raw material (sodium bicarbonate). The powder from the *bibag* meets the requirements of the US Pharmacopoeia. The safety and efficacy of the *bibag* product are not affected by these discolorations.

STORAGE: Store between 5°C and 30°C (41°F and 86°F). Do not refrigerate. In order to avoid any contamination, keep the *bibag* tightly sealed until it is connected to the hemodialysis machine.

CHEMICAL COMPOSITION OF CONCENTRATE: Sodium Bicarbonate, USP.

PRESENTATION: 650g or 900g *bibag* bicarbonate concentrate in polyamide/polyethylene container.

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