

# FRESENIUS MEDICAL CARE

## Kinexus<sup>®</sup> Wi-Fi Gateway for Liberty User's Guide

For use with the Liberty<sup>®</sup> Select Cyclor

P/N 480173 Rev D | July 2024

## Legal

© 2024 Fresenius Medical Care. All Rights Reserved.

Fresenius Medical Care, the triangle logo, Liberty, Liberty Select, and Kinexus are trademarks of Fresenius Medical Care Holdings, Inc. or its affiliated companies. All other trademarks are the property of their respective owners.

This document contains proprietary information of Fresenius Medical Care Renal Therapies Group, LLC and its affiliates (“Fresenius Medical Care”). The contents of this document may not be disclosed to third parties, copied, or duplicated in any form, in whole or in part, without the prior written permission of Fresenius Medical Care.

Fresenius Medical Care only distributes the Gateway for use with the Liberty Select Cyclor. For more information, contact:

Fresenius USA, Inc.  
Waltham, MA 02451  
(800) 227-2572

This document provides overview, installation instructions, specifications, and troubleshooting for the Gateway. For more information about the Liberty Select Cyclor, refer to the *Liberty Select 3.0 User's Guide*, P/N 480165.

**Use:** The Kinexus Portal is intended to allow clinicians to review patient data collected from the patient’s home dialysis device(s) and to enable physicians to manage patients’ home dialysis therapies for certain Fresenius Medical Care products, including peritoneal dialysis cyclor(s). The Kinexus Portal does not provide medical diagnosis or recommendations regarding medical treatment. The Kinexus Portal and the content appearing on this site are not substitutes for professional medical advice, diagnosis, or treatment.

The Kinexus Gateway is intended for transmission of home dialysis therapy data and patient data for patients using certain Fresenius Medical Care products, including peritoneal dialysis cyclor(s).

Kinexus is not an electronic health record (EHR) or billing solution.

**Indications for Use:** The Fresenius Medical Care Liberty Select Cyclor is indicated for acute and chronic peritoneal dialysis.

Patents apply, see [www.fmcna.com/patents](http://www.fmcna.com/patents)

# TOC - 1

---

## Glossary 1

### Chapter 1: Introduction

<b>Introduction</b> .....	<b>1-1</b>
About this Guide - Call-outs .....	1-1
Using this Guide .....	1-2
What's Included .....	1-4

### Chapter 2: Gateway Overview

<b>Overview</b> .....	<b>2-1</b>
Front View .....	2-2
Back View .....	2-3
Top View .....	2-4
Bottom View .....	2-6
Power Supply .....	2-7

### Chapter 3: Gateway Setup

<b>Gateway Setup</b> .....	<b>3-1</b>
<b>Connect Gateway to Your Home Wi-Fi</b> .....	<b>3-4</b>
Internet .....	3-4
Installing and Launching the Kinexus Wi-Fi App for Android Users .....	3-5
Replacement Scenario For Android Users .....	3-19
Installing and Launching the Kinexus Wi-Fi App for iOS Users .....	3-28
Replacement Scenario for iOS Users .....	3-39
Connect Gateway to Cyclor .....	3-47
Gateway Commands .....	3-51

### Chapter 4: Troubleshooting

Gateway .....	4-1
Android Application .....	4-2
iOS Application .....	4-4

## TOC - 2

---

### Chapter 5: Appendix

<b><i>Device Specifications</i></b> .....	<b>5-1</b>
Dimensions.....	5-1
Operating Conditions.....	5-1
Storage Conditions.....	5-1
Maximum Altitude.....	5-1
Electrical.....	5-2
Regulatory.....	5-3
FCC Information.....	5-3
WEE Information.....	5-4
<b><i>Radio Bands and Frequencies</i></b> .....	<b>5-5</b>
North America Modem & Wireless.....	5-5
Bluetooth.....	5-5
Wi-Fi.....	5-5
Marker Plate Symbols.....	5-6
Cleaning Your Gateway.....	5-7

### Chapter 6: Gateway Information

<b><i>Gateway Information</i></b> .....	<b>6-1</b>
---	------------



## Glossary

---

**Home Router:** Home router will be connected by the Kinexus mobile app to the Wi-Fi Gateway.

**Cloud:** A place on the internet where your Gateway and your healthcare team can exchange information. Data on the cloud is kept secure.

**Ethernet Cable:** A type of cable that connects two devices through a local area network (LAN).

**LAN:** A type of port connected by an Ethernet cable. Ethernet cables connect devices to each other or connect a device to the internet. Stands for *Local Area Network*.

**LED:** A small light on the top of the Gateway. Stands for *Light Emitting Diode*.

**RS-232:** A type of port used to send data between two devices. The RS-232 port on the Liberty Select Cyclor connects to a USB port on the Gateway. Stands for *Recommended Standard 232*.

**USB:** A type of port and cord. USB cords connect devices to each other. Stands for *Universal Serial Bus*.

**END OF CHAPTER**

# Chapter 1: Introduction

## Introduction

Communication between you and your health care team is essential to making sure you get the most out of your dialysis treatments. The Kinexus Gateway for the Liberty Select Cyclor (“Gateway”) is a new accessory made to assist your healthcare team in keeping track of your Peritoneal Dialysis (PD) therapy. The Gateway transfers treatment and alarm history data at the end of each treatment to your care team. A home router is required to use the Gateway.

The Gateway may receive therapy parameters from your healthcare team. The Gateway is not intended to upload any software or treatments to the Liberty Select Cyclor and is not a treatment monitoring device.

## About this Guide - Call-outs

You should read this guide before using the Gateway with your cyclor. The guide has several helpful call-outs within to help you use the Gateway correctly.



### WARNING

Warnings warn against actions that could hurt you. If you do not follow these directions, it can lead to serious injury or death.



### WARNING - SHOCK HAZARD

Shock Hazard means that you could get a severe electrical shock if you do not use the equipment the right way.

**CAUTION**

Cautions warn against actions that could harm you or your machine.

---

**NOTE:**

Notes mark extra information about Gateway features. Notes have helpful information. They may also have information about unique cases.

---

## Using this Guide

This user's guide is designed to walk you through the process of setting up your Gateway to work with your cyclor. Look through the Table of Contents to quickly find a section, or search the Index to find a specific topic. Pages are also labeled on the sides to help you keep track of what section you are reading.

**NOTE:**

The Gateway is not in use during the treatment. Some wireless transmissions from the Gateway may be actively emitted during treatment or device functions. The associated hazards should be minimal.

---

**NOTE:**

If the treatment data transfer is not successful, your cyclor will retry the transfer at the end of next treatment.

---

**NOTE:**

Do not connect any other accessories to your cyclor that are not approved by Fresenius Medical Care.

---

**NOTE:**

Do not modify the Gateway provided to you by Fresenius Medical Care.

---

**NOTE:**

Do not connect or disconnect the Gateway from your cyclor when the cyclor is powered ON.

---



## What's Included

When you receive your Gateway, you will find the following items:

- The Gateway
- One (1) 100-240V AC wall power adapter
- One (1) RS-232 to USB cable
- Kinexus Gateway User's Guide (this document)
- SysLINK™ / SysCoRE™ Gateway User Manual and Installation Guide

**NOTE:**

Images may differ from actual items.

---

**NOTE:**

For all inquiries and questions, call Fresenius Medical Care Technical Support at **(800) 227-2572**.

---

**END OF CHAPTER**



## Chapter 2: Gateway Overview

---

### Overview

The Gateway is a small device that connects to the internet to communicate with your healthcare team through the cloud. The Gateway communicates with the cloud to transfer treatment and alarm history data at the end of each treatment. The Gateway has ports to connect the device. There are several LED lights on top of the Gateway. These lights help show you what the device is doing or if there is a problem with the Gateway.



#### **CAUTION**

---

Do not share your Gateway or give your Gateway to others to use. The Gateway must be distributed to and used only by the person it is prescribed for.

---

## Front View




*Front View*

The front of the Gateway is a blue plate with the Fresenius Medical Care logo on the left. There are no LEDs or ports on the front.

The left side of the Gateway is a heat vent. The right side is a speaker. The Gateway lets you know about its status with voice responses. See [Gateway Commands](#) on [page 3-51](#) for a list of the possible responses you may hear.

## Back View



*Back View*

The back of the Gateway is where you power on the Gateway and where you connect the Gateway to your cyclor. The back of the Gateway includes:

- 1 Power Cord Inlet:** This is where you plug in the power cord. Plugging in the power cord turns the Gateway ON. Unplugging the power cord turns the Gateway OFF.
- 2 Reset button:** The Reset button is located on the back panel of the gateway. To reset the gateway, press the Reset button momentarily (less than one second). Use a paper clip, or similar object, to press. The Status light will turn orange initially during the reset.

To reset the Gateway back to factory defaults, press and hold the Reset button until the status light alternates between red and green. This will occur in about 10 seconds. Then release the Reset button.




### CAUTION

Do not press reset your Gateway without guidance from Technical Services support. Pressing this button will reset your Gateway and your Gateway will no longer function.

- 3 USB 1 & 2 port:** This is where you plug in the USB-to-RS232 cable that connects the Gateway to the cyclor. Kinexus Gateways have either (1) or (2) USB ports. Either port may be used to connect your Liberty Select cyclor.
- 4 WAN/LAN Port:** The WAN/LAN Port is not used.
- 5 LAN Port:** The LAN Port is not used.

## Top View

The top of the Gateway displays several LEDs. These LEDs show the power status, if the Gateway is communicating with the cyclor.



*Top View*

The top parts of the Gateway, include:

- 1 USB 1 & 2:** The USB1 and USB2 LED lights are green when connected; red when a fault has occurred, amber/orange when disabled, and off when not in use.
- 2 Option:** No function.
- 3 Bluetooth:** The *Bluetooth* LED will be solid green or flicker green during normal operation. The light is off if the Bluetooth is not connected or not functioning. If DISABLED (amber/ orange color), the feature must be enabled within the *Setup* Menu.
- 4 Wi-Fi:** The LED is solid green when Wi-Fi is connected; when Wi-Fi is not connected, this LED is orange.
- 5 WAN/LAN:** No function.

- 6 Status:** The *Status* LED indicates current gateway status. It will be solid amber/orange at initial power up. The light will blink green at one second intervals during normal operation. Should the status light stop blinking, remaining in the solid OFF or ON state for more than 60 seconds, remove power supply from the wall plate power outlet, wait 30 seconds and then restore power by plugging the power supply into the wall plate power outlet. Wait until the status light begins blinking again before using the Gateway.

These LEDs show the status of the connection.

- Solid green means that the Gateway is ON.
- Flashing green means that the Gateway is receiving or transmitting data.
- Flashing red and green means that the Gateway is not connected to the internet.

- 7 Setup:** No function.

## Bottom View



*Bottom View*

The bottom panel of the Gateway is labelled with the marker plate. The marker plate lists the manufacturer name, model number, part number, serial number, and any cautions or warnings associated with the device. Some of the information on the bottom of the Gateway is unique to each Gateway and your Gateway will not exactly match the image above.

For more information, see [Marker Plate Symbols](#) on [page 5-6](#).

## Power Supply



Plug the round power supply connector to the Gateway. Then plug the power supply directly into a wall plate power outlet.



### **WARNING - SHOCK HAZARD**

To power on the gateway, always insert the round power supply connector to the gateway first. Then insert or plug in the power supply to the wall plate power outlet.

---



### **WARNING - SHOCK HAZARD**

To power off the Gateway, always remove the power supply from the wall plate power outlet first. Then remove the round power supply connector from the Gateway.

---

**END OF CHAPTER**



## Chapter 3: Gateway Setup

### Gateway Setup

Before setting up your Gateway, you should find a flat, dry surface near your cyclor to set it up. Make sure the Gateway is at least eight (8) inches away from your body during your dialysis treatment.

**NOTE:**

Avoid placing the Gateway near your bedside if at all possible. Electromagnetic (EM) radiation may interfere with your other medical devices.

---

**CAUTION**

Pacemaker devices, implanted neurostimulators, and implanted defibrillators may be affected by EM radiation. If you wear a pacemaker, implanted neurostimulator, or implanted defibrillator, you should avoid contact with the Gateway and stay further away from the Gateway than the recommended eight inches. You should also avoid installing the Gateway yourself.

---

**CAUTION**

Position the power cord and RS232 to USB cable so they will not be tripped over by you or others living in your home.

---

**CAUTION**

Do not let the Gateway come into contact with blood or bodily fluids, i.e., dialysate from your abdomen.

---

**CAUTION**

Place the Gateway on a flat, dry surface, away from any potential fluid spills or splatter.

---

**NOTE:**

The Gateway may transmit during treatment or other cyclor functions.

---

**NOTE:**

The wall plate power outlet should be near to the Gateway and easily accessible.

---

**NOTE:**

Physical placement of the Gateway near a bedside should be avoided where possible due to EMC radiation.

---

**CAUTION**

Pacemaker devices may be affected by EMC in several ways and wearers should avoid contact with the Gateway. All persons should remain at least 20 cm (8 inches) from the transmitter's antenna. Pacemaker wearers should try to increase this distance where possible for general safety margins of operation and avoid self-installation of the Gateway.

---

Your Gateway may communicate with other devices using wireless radio, Ethernet, or a serial/USB cable. Your device supplier will furnish you with installation instructions specific to your devices.

Generally, your additionally supplied devices will be either for medical or information technology applications. In some cases, the two applications can be combined. When that occurs, follow the medical application instructions where applicable.



**WARNING - SHOCK HAZARD**


If your Gateway is being used around medical patients, the Gateway must be out of reach by the patient. Some patients may have a weakened immune system.

---


# Connect Gateway to Your Home Wi-Fi

## Internet

- 1 Remove your Gateway from the packaging. Inspect the Gateway for signs of damage such as breaks or cracks.
- 2 Plug the round end of the power cord into the power cord inlet in the back of the Gateway, **2**.




- 3 Plug the power adapter directly into a 120v wall outlet.
- 4 Kinexus Gateways have either (1) or (2) USB ports. Either port may be used to connect your Liberty Select cyclor.
- 5 Once the Gateway has been plugged in the Status LED will turn orange. Wait until the LED turns green to continue.




## Installing and Launching the Kinexus Wi-Fi App for Android Users

**Prerequisite:** The Wi-Fi Gateway requires a smart phone or tablet to connect to your home internet service. If your health care provider issues you a mobile device for use with your therapy systems this application may already be installed. Please use your issued mobile device to setup the gateway if the application is already present.


- 1 Enter “**Kinexus Wi-Fi**” in the search bar of the Google Play Store.
- 2 Select the Kinexus Wi-Fi application.
- 3 Click on the “**Install**” button to install the Kinexus Wi-Fi app.
- 4 After the successful installation, you will be presented with a screen similar to the one shown in the figure below. To launch the application, click on the “**Open**” button.




5 When prompted for location access, select the “While using the app” option.



- 6 Select “**Allow**” when prompted to find nearby devices, as shown in the figure below.




- 7 Click on the “**START**” button to initiate the search for nearby Gateways.







- 8 The Kinexus app will begin searching for nearby Gateways, as shown in the figure below.




- 9 Once the application has identified the available Gateways, you can choose your preferred Gateway from the list provided below.




- 10 After selecting your preferred Gateway from the list, the application will begin the process of connecting to the Gateway, as shown in the figure below.




- 11 Once the Gateway is successfully connected, the application will initiate a scan to detect available Wi-Fi networks. Please note that this scanning process may take several minutes to complete.




- 12 A list of available Wi-Fi networks will be displayed, as shown in the figure below. From this list, select your home Wi-Fi network.



13 After selecting your home Wi-Fi network, click on the “**CONTINUE**” button.




- 14 Enter the password of your home Wi-Fi network and click on the “**APPLY**” button.




The screenshot shows a mobile application interface for Kinexus. At the top, the Kinexus logo is displayed. Below the logo, the text reads "Please enter password for Abracadabra6-5:". There is a white rectangular input field for the password. Below the input field, there is a checkbox labeled "show password". At the bottom of the screen, there are three buttons: "APPLY", "BACK", and "START OVER". The bottom of the screen shows the standard Android navigation bar with three icons: a home button, a back button, and a recent apps button.

- 15 Once the password is entered, the app will initiate the connection to the Wi-Fi network. Please note that the process of connecting to the network may take several minutes to complete.






- 16 If the entered password is incorrect or if the Gateway fails to establish an internet connection, you will see the screen displayed in the figure below. In this case, please check your internet connection, re-enter your password, and click on the “**APPLY**” button.




- 17 Once the setup of your Gateway is completed and it successfully connects to your home Wi-Fi network, you will see the Success screen displayed, as shown in the figure below. This indicates that the Gateway setup process has finished successfully.




## Replacement Scenario for Android Users


- 1 After the application has identified the available Gateways, you have the option to select a different Gateway if needed. From the list provided below, you can choose your preferred Gateway by selecting it.




- 2 After selecting your preferred Gateway from the list, the application will begin the process of connecting to the Gateway, as shown in the figure below.




- 3 Once the Gateway is successfully connected, the application will initiate a scan to detect available Wi-Fi networks. Please note that this scanning process may take several minutes to complete.




- 4 A list of available Wi-Fi networks will be displayed, as shown in the figure below. From this list, select your home Wi-Fi network.



- 5 After selecting your home Wi-Fi network, click on the “**CONTINUE**” button.




- 6 Enter the password of your home Wi-Fi network and click on the “**APPLY**” button.







- 7 Once the password is entered, the app will initiate the connection to the Wi-Fi network. Please note that the process of connecting to the network may take several minutes to complete.



- 8 If the entered password is incorrect or if the Gateway fails to establish an internet connection, you will see the screen displayed in the figure below. In this case, please check your internet connection, re-enter your password, and click on the “**APPLY**” button.




- 9 Once the setup of your Gateway is completed and it successfully connects to your home Wi-Fi network, you will see the Success screen displayed, as shown in the figure below. This indicates that the Gateway setup process has finished successfully.




## Installing and Launching the Kinexus Wi-Fi App for iOS Users

**Prerequisite:** The Wi-Fi Gateway requires a smart phone or tablet to connect to your home internet service. If your health care provider issues you a mobile device for use with your therapy systems this application may already be installed. Please use your issued mobile device to setup the gateway if the application is already present.


- 1 Enter “**Kinexus Wi-Fi**” in the search bar of the Apple App Store.
- 2 Select the Kinexus Wi-Fi application.
- 3 Click on the “**Install**” button to install the Kinexus Wi-Fi app.
- 4 After the successful installation, you will be presented with a screen similar to the one shown in the figure below. To launch the application, click on the “**Open**” button.




- 5 Click on the “**START**” button to initiate the search for nearby Gateways.




- 6 The Kinexus app will begin searching for nearby Gateways, as shown in the figure below.



- 7 Once the application has identified the available Gateways, you can choose your preferred Gateway from the list provided below.




- 8 After selecting your preferred Gateway from the list, the application will begin the process of connecting to the Gateway, as shown in the figure below.







- 9 When prompted with a Bluetooth Pairing Request, select the “**Pair**” option as shown in the figure below.




- 10 Once the Gateway is successfully connected, the application will initiate a scan to detect available Wi-Fi networks. Please note that this scanning process may take several minutes to complete.



- 11 A list of available Wi-Fi networks will be displayed, as shown in the figure below. From this list, select your home Wi-Fi network.




- 12 Enter the password of your home Wi-Fi network and click on the “**APPLY**” button.




The screenshot shows a mobile application interface for Kinexus. At the top, the time is 19:42 and there are icons for cellular signal, Wi-Fi, and battery. The Kinexus logo is displayed below the status bar. The main content area has a light gray background and contains the text "Please enter password for Abracadabra6-2:". Below this text is a white password input field. Underneath the input field is a checkbox labeled "show password". At the bottom of the form are three buttons: "APPLY", "BACK", and "START OVER".

- 13 Once the password is entered, the app will initiate the connection to the Wi-Fi network. Please note that the process of connecting to the network may take several minutes to complete.




- 14 Once the setup of your Gateway is completed and it successfully connects to your home Wi-Fi network, you will see the Success screen displayed, as shown in the figure below. This indicates that the Gateway setup process has finished successfully.




## Replacement Scenario for iOS Users

- 1 After the application has identified the available Gateways, you have the option to select a different Gateway if needed. From the list provided below, you can choose your preferred Gateway by selecting it.




- 2 After selecting your preferred Gateway from the list, the application will begin the process of connecting to the Gateway, as shown in the figure below.







- 3 When prompted with a Bluetooth Pairing Request, select the “**Pair**” option as shown in the figure below.




- 4 Once the Gateway is successfully connected, the application will initiate a scan to detect available Wi-Fi networks. Please note that this scanning process may take several minutes to complete.



- 5 A list of available Wi-Fi networks will be displayed, as shown in the figure below. From this list, select your home Wi-Fi network.




- 6 Enter the password of your home Wi-Fi network and click on the “**APPLY**” button.



- 7 Once the password is entered, the app will initiate the connection to the Wi-Fi network. Please note that the process of connecting to the network may take several minutes to complete.



- 8 Once the setup of your Gateway is completed and it successfully connects to your home Wi-Fi network, you will see the Success screen displayed, as shown in the figure below. This indicates that the Gateway setup process has finished successfully.




## Connect Gateway to Cycler


**NOTE:**

Do not connect or disconnect the Gateway from your cycler while the cycler is powered ON.


- 1 With the Power Switch on your Liberty Select Cyclor set to OFF, connect the USB-to-RS-232 cable to your Liberty Select Cyclor, **a**.
  - a. If your Liberty Select Cyclor has a cap over the RS-232 port, remove the RS-232 cap using a screwdriver.



- b. If your Liberty Select Cyclor does not have a cap over the RS-232 port, proceed to step 1c.






c. Connect the RS-232 end of the RS-232 to USB cable to the RS-232 port on your Liberty Select Cyclor, .



- 2 Plug the USB end of the RS-232 to USB cable into the USB port on the back of the Gateway, **g2**. Kinexus Gateways have either (1) or (2) USB ports. Either port may be used to connect your Liberty Select cyclor.



## Gateway Commands

**Failed to upload your treatment data** – The Gateway was unable to upload your treatment information to the cloud. It will re-try automatically at the end of your next treatment.

**Gateway is ON** – Confirms the Gateway is powered ON.

**Sending your treatment information to your clinic, please wait** – Confirms the Gateway successfully received your treatment information from your cyclor and is transmitting data to the cloud.

**There seems to be a connectivity issue with your Gateway. Please contact technical support using the number found on your cyclor** – Your Gateway does not have an internet connection.

**Treatment data uploaded successfully** – Confirms the Gateway successfully uploaded your treatment information to the cloud.

END OF CHAPTER

## Chapter 4: Troubleshooting


### Gateway

**NOTE:**

The Gateway is non-user-serviceable and must be replaced by Fresenius Medical Care. Call Fresenius Medical Care Technical Support at **(800) 227-2572** for questions regarding set-up, troubleshooting, or assistance in returning or replacing your Gateway.

If no LEDs light up on your Gateway:


- 1 Unplug the power cord from the Gateway and the wall outlet. Plug the power cord back into the Gateway, and then plug the cord back into the wall outlet.
- 2 If the outlet is connected to a power switch, ensure that the switch is set to ON.
- 3 Check to make sure the power to the wall outlet is working.
  - a. Unplug the Gateway from the wall outlet and plug another device, like a lamp, into the same outlet.
    - If device works, unplug it and plug the Gateway back into the wall outlet.
    - If the device doesn't work, use a different outlet to test the Gateway.
- 4 If the outlet is a GFI (ground-fault interrupter):
  - a. Check to see if the Reset button is popped out on the outlet. Reset buttons are usually the large button in the middle of the wall outlet, and are often colored red. See the diagram on the right for an example of a GFI outlet.
  - b. If the Reset button is popped out, unplug the Gateway from the wall outlet, press the Reset button back in, and plug the Gateway back into the outlet.
- 5 Check your fuse box to see if a circuit breaker needs to be reset.



If none of these steps work, call Fresenius Medical Care Technical Support at **(800) 227-2572**.


## Android Application

If you encounter any issues during the setup of your Gateway using the Kinexus Wi-Fi application on your Android mobile device, please contact Fresenius Medical Care Technical Support at **+1 (800) 227-2572** for assistance.



When you encounter the screen shown below upon launching the Kinexus Wi-Fi application on your Android mobile device, following any of the steps mentioned below can help resolve the issue:


- 1 Check that the Gateway is properly plugged in and connected.
- 2 Ensure that your mobile device is in close proximity to the home router for a strong and stable signal.
- 3 Check the Wi-Fi LED status, and if it is solid green, then Gateway is already connected to your home Internet service. In this case, you can exit the app.




## iOS Application

When you encounter the screen shown below upon launching the Kinexus Wi-Fi application on your iOS mobile device, following any of the steps mentioned below can help resolve the issue:

- 1 Check that the Gateway is properly plugged in and connected.
- 2 Ensure that your mobile device is in close proximity to the Gateway for a strong and stable connection.
- 3 Check the Wi-Fi LED status, and if it is solid green, then Gateway is already connected to your home Internet service. In this case, you can exit the app.




When you encounter the screen shown below upon launching the Kinexus Wi-Fi application, please ensure that Bluetooth is enabled on your iOS mobile device.






If you encounter any issues during the setup of your Gateway using the Kinexus Wi-Fi application on your iOS mobile device, please contact Fresenius Medical Care Technical Support at **+1 (800) 227-2572** for assistance.



If the entered password is incorrect or if the Gateway fails to establish an internet connection, you will see the screen displayed in the figure below. In this case, please check your internet connection, re-enter your password, and click on the **“APPLY”** button.



**END OF CHAPTER**



## Chapter 5: Appendix

---

# Device Specifications

### Dimensions

**Floor space:** 5.5 inches (13.97 cm) wide by 4 inches (10.16 cm) deep

**Height:** 1 inch (2.54 cm)

**Weight:** 1lb/0.45kg

### Operating Conditions

**Temperatures:** 50° — 104° F (10° C — 40° C) with provided power supply. For use at environments below 10°C or above 60°C, use a power source rated for 60°C, 9-24VDC, minimum 12W, tip positive, listed to UL/ANSI 60601-1 and UL/ANSI 60950-1.

**Relative Humidity:** 10%-85% (non-condensing)

### Storage Conditions

**Temperatures:** -22° — 149° F (-30° C — 70° C)

**Relative Humidity:** 5%-90% (non-condensing)

### Maximum Altitude

9842ft/3000m

## Electrical

**Power Supply:** Positive polarity Direct Current at 9 - 24V and minimum 18W input capacity connected to a circuit equipped with a hospital grade receptacle

**External Connections:** Ethernet ports; USB port

## Regulatory

Your Gateway was tested for operation in an environment such that medical devices and other medical equipment can be used. Your Gateway meets IEC 60601-1-2 and FCC Class B requirements for use in business, industrial, commercial, and residential environments. In addition, your Gateway meets the safety standards for IEC 62368-1.

The Marker Plate located on the bottom of the Gateway will contain regulatory information about radios installed.

The Gateway complies to:

Colocation – FCC 15.247 (FCC Part 15 Subpart C)

Transmitters Colocation Evaluation Applicable to

USA & Canada

EU

Australia/New Zealand

FCC15B (FCC Part 15 Subpart B)

IEC 62368-1, second and third edition

IEC 60601-1-2, 4th edition

## FCC Information

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Increase the separation between the equipment and cyclor.
- Connect the equipment into an outlet on a circuit different from that to which the cyclor is connected.

## WEE Information

This is important information for users with regards to the proper disposal and recycling of your Gateway. You are required to comply with this notice for all electronic products bearing the following symbol:



# Radio Bands and Frequencies

## North America Modem & Wireless

Module Quectel SC20-AL / WL

Modem

Frequencies:

### Bluetooth

BT 2.1+EDR/3.0/4.1 LE/4.2 BLE

2.4GHz (2402-2480 MHz)

### Wi-Fi

802.11a/b/g/n

2.4Ghz (2400-2482 MHz)

5Ghz (5180-5825 MHz)

## Marker Plate Symbols

Your Marker Plate is the label located on the bottom of your Gateway. The Marker Plate contains important information as to the environmental conditions for usage. The symbols below define the conditions of usage required for the gateway and the patient.



This symbol indicates that accompanying documentation with the Gateway must be read prior to installation.



This symbol indicates the temperature range.



This symbol indicates the humidity range.




This symbol indicates the elevation range.



This symbol indicates the weight of the Gateway.



This symbol indicates the power input to the Gateway is DC power at the defined voltage and a minimum watt input capability.



This symbol indicates positive polarity for the DC power input to the Gateway.



## Cleaning Your Gateway

The Gateway is not sterile equipment; however, you should keep your Gateway clean and dry at all times. Clean the **exterior surface** of your Gateway using a damp cloth and a diluted (1:100) bleach solution or other suitable hospital disinfectant.

Use surface cleaner in **small amounts** to avoid excess cleaner from entering the interior of your Gateway. Rinse off cleaning solutions with a water-dampened cloth, especially if a corrosive cleaner such as bleach is used.

---



### CAUTION

---

Do not allow any fluid to drip into your Gateway. The Gateway must be OFF when you are cleaning it. If any fluid gets inside your Gateway, you must contact a Fresenius Medical Care Technical Support immediately at **(800) 227-2572**. Clean fluid spills immediately to prevent moisture from possibly damaging your Gateway. Your Gateway should be placed on a water-tolerant surface, and spills should be cleaned up immediately.

---



### CAUTION

---

Do not use foaming type cleansers or disinfectants containing alcohol or quaternary ammonium compounds like N-alkyl (C12-C18) dimethyl benzyl ammonium chloride (benzalkonium chloride). These chemicals attack the polycarbonate plastics used in your machine. Read the product labels and follow the instructions. Do not use spray bottles.

---

Any service must be performed by a Fresenius Medical Care technician. If your Gateway needs service, call Fresenius Medical Care Technical Support at **(800) 227-2572**.

**END OF CHAPTER**



## Chapter 6: Gateway Information

---

### Gateway Information

Please complete the following information for future reference. You will find most of this information on the bottom label of your gateway and your power supply.

*Serial Number:* \_ \_ \_ \_ \_

*Model Number:* SL-600 IP22 / SC-600 IP 22

*MEID:* \_\_\_\_\_

*Internet Connection:*  Wi-Fi

END OF CHAPTER

---

This page intentionally left blank

---

This page intentionally left blank