PREVENTIVE MAINTENANCE CHECKLIST ANNUAL/4000 HOUR

SERIAL NO.		DATE			
MAC					
TECH	HNICIAN(S)				
	See the Referenced Section of the 2008K@home [™] Preventive I (part number 507665 Rev. K) for detailed instructions to perfor			S	
Refere	ence Section Procedure		Value	Completed By:	
3.1	Clean Filters and Replace O-Rings				
3.2	Replace bibag Inlet Filter and Connector Door O-Rings				
3.3	Replace Pre-UF Pump Filter				
3.4	Replace UF Pump Check Valves				
3.5	Check Deaeration Restrictor and Replace O-Rings				
3.6	Rebuild UF, Concentrate and Bicarbonate Pumps				
3.7	Check Ground Heater Element				
3.8	Check High Voltage AC Connections			-	
3.9	Replace Deaeration Motor Brushes(e				
3.10	Check Inlet Water Pressure Regulator				
3.11	Calibrate bibag Pressure Transducers				
3.11	Check Deaeration Pressure				
3.12	Check Loading Pressure (Range: 18 to 20psi / with DIASAFE* Plus				
3.13					
	Check Flow Relief Pressure (Range: 29 to 30psi/with DIASAFE* Plus				
3.14	Test Concentrate and Bicarbonate Pumps				
3.15	Calibrate UF Pump (Range: 24 ± 0.1cc at 24 U	F Pump strokes) _		· -	
3.16	Check Conductivity Calibration(Range: Δ ± 0.1 mS/cm)		External Meter		
3.17	Check Temperature Calibration(Re	ange: 37 ± 0.3°C)			
	(Ra	ange: 39 ± 0.3°C)			
3.18	Check Volt Hi Lo Detect(Range:				
3.19	Check Blood Leak Calibration(R				
	Check Dimness Calibration(R	ange 4.0 to 6.0V)			
3.20	Check Arterial, Venous and TMP Calibration	•••••			
3.21	Check Dialysate Flow(Range @ 500ml/m	in: 485 to 515ml)			
3.21	Check Heparin Pump	•• (If Applicable)			
3.22	Check Blood Pump	•••••			
3.23	Check and Calibrate Level Detector	•••••			
3.24	Check Alarm Operation and Pressure Holding Tests	•••••			
3.25	Verify pH	(Range 6.9 to 7.6)			
3.26	Rinse Check(Range: >	310ml collected)			
	Heat Disinfect Check(Ra				
3.27	Check Power Failure Alarm	_			
	Replace 9-Volt Battery	•••••			
3.28	Check Blood Pressure Module				
3.29	Final Checks				