

2008T BlueStar™ Hemodialysis Machine

Training Guide



**FRESENIUS
MEDICAL CARE**

RENAL TECHNOLOGIES

Indication for Use:

2008T BlueStar Hemodialysis Machine: The 2008T BlueStar Hemodialysis Machine is indicated for acute and chronic dialysis therapy in a healthcare facility.

Additional therapy options for patients receiving hemodialysis include: Isolated Ultrafiltration, Sustained Low Efficiency Dialysis (SLED), and low volume hemodialysis (patients weighing ≥ 20 kg and ≤ 40 kg). This machine accommodates the use of both low flux and high flux dialyzers. The SLED therapy option is not to be used for patients weighing ≤ 40 kg. The 2008T BlueStar Hemodialysis Machine is not to be used for plasma replacement therapies, for patients weighing less than 20 kg, or for renal therapies using substitution fluid.

Caution: Federal (US) law restricts this device to sale by or on the order of a physician.

Note: Read the Instructions for Use for safe and proper use of this device. For a complete description of hazards, contraindications, side effects and precautions, see full package labeling at fmcna.com.

2008T BlueStar Hemodialysis Machine Training Guide

The *2008T BlueStar™ Hemodialysis Machine Training Guide* is designed for the training of skilled hemodialysis nurses and technical assistants in the use of the machine. Use this guide as a companion to the *2008T Hemodialysis Machine Operator’s Manual*, which contains detailed instructions for all machine functions, and the *2008T Hemodialysis Machine Quick Start Guide*.

This Training Guide covers the basic instruction topics for the 2008T BlueStar Hemodialysis Machine. It is not intended to compromise or replace specific policies and procedures of individual hemodialysis providers. Fresenius Renal Technologies and its affiliates make no recommendations and assume no liability for the policies and procedures established in hemodialysis provider’s clinics.

The guide is designed to be flexible. While approximate timeframes for training sessions are given, the actual time commitment required will depend on the size of the group and the prior experience of participants.






Similarly, the guide suggests topics where hands-on exercises will enhance your participant’s experience. These demonstrations may occur at many different times during your facility’s training, with the time and location dependent on the individual unit. We encourage you to tailor these hands-on experiences to the needs of the group and to pace the sessions accordingly.

In addition to this Training Guide, the following materials are integral to training in the proper use of the 2008T BlueStar Hemodialysis Machine:

- *2008T Hemodialysis Machine Operator’s Manual (P/N 490122)*
- *2008T Hemodialysis Machine Quick Start Guide (P/N 490161)*
- *2008T Operator’s Troubleshooting Guide with bibag (P/N 102297-01)*
- *2008T BlueStar Hemodialysis Machine Participant Workbook (P/N 103324-01)*

Icons and Descriptions:

The following icons are used within this Training Guide to provide direction and clarity:

Icons and Descriptions	
	Trainer explanation and demonstration
	Facilitation tip and supplemental comments
	Participant hands-on practice recommended
	Information found in the <i>Quick Start Guide</i>
	Features unique to the 2008T BlueStar Hemodialysis Machine

Training Timeline

3-4 weeks prior to training:

- Determine the number of participants.
- Order training supplies to be shipped directly to facility (2008T BlueStar Machine Getting Started Training Kit).

2 weeks prior to training:

- Confirm dialysis supplies are on hand or have been ordered.

1 week prior to training:

- Ensure machines have arrived and are installed.
- Ensure cultures have been drawn.
- Ensure dialysis supplies are in the unit.
- Confirm receipt of the training materials.

2008T BlueStar Machine Getting Started Training Kit, includes:

- *2008T Hemodialysis Machine Operator's Manual (P/N 490122)*
- *2008T BlueStar Hemodialysis Machine Training Guide (P/N 103323-01)*
- *2008T Hemodialysis Machine Quick Start Guide (P/N 490161), (one per participant)*
- *2008T BlueStar Hemodialysis Machine Participant Workbook (P/N 103324-01), (one per participant)*
- *2008T Operator's Troubleshooting Guide (P/N 102297-01)*
- *2008T BlueStar In-Service Checklist (P/N 103322-01)*

Training Topics and Agenda

	Estimated training time (minutes)
Welcome/Getting Started.....	5
<ul style="list-style-type: none"> • <i>Participant Workbook</i> Distribution and Introduction • <i>2008T Hemodialysis Machine Operator's Manual</i> Introduction • <i>2008T Hemodialysis Machine Quick Start Guide</i> Introduction 	
Review 2008T BlueStar Hemodialysis Machine.....	10-15
<ul style="list-style-type: none"> • Review Key Machine Features 	
Demonstrate Daily Preparation for Treatment.....	45
<ul style="list-style-type: none"> • Set Up Concentrate • Set Up the Bloodline • Prime the Extracorporeal Circuit, Manual, or Auto Prime • Test the 2008T BlueStar Hemodialysis Machine • Recirculation • Manual or Auto Prime • Manual or Auto Start 	
Review Screens	
Review Initiation of Treatment.....	10
Alarms and Troubleshooting Activity.....	15
Review Termination of Treatment.....	5-10
<ul style="list-style-type: none"> • Manual or Assisted Reinfusion 	
Disinfection and Maintenance.....	10-15
<ul style="list-style-type: none"> • Acid Clean and Heat Disinfect • Acid Clean and Chemical Rinse • DIASAFE® Plus Filter Replacement • Manual or Assisted Reinfusion 	
Appendix	
<ul style="list-style-type: none"> • HD Equipment and MIS Support • In-Service Evaluation 	

Getting Started



Welcome Participants

Have all of your registered participants sign the In-Service Attendance Sheet. Solicit the cooperation of the group to ensure a good learning experience.

Set expectations:

- No cell phone usage, including text messaging.
- Show respect for other participants.
- Share all comments and questions with the room. (i.e., no side conversations)
- Consider using a “parking lot” for questions, particularly if the group is large.

Distribute and introduce the *2008T BlueStar Hemodialysis Machine Participant Workbook* (P/N 103324-01)

- Give participants a few minutes to familiarize themselves with the workbook. Encourage participants to take notes for future reference after the training.
- Introduce the In-Service Checklist, on page 4 of the Participant Workbook, and encourage participants to check topics off as they are covered during the training. Individual facilities may want to keep a copy of the Checklist in the employee’s file.

Introduce the *2008T Hemodialysis Machine Operator’s Manual* (P/N 490122)

- Point out the Alarms/Troubleshooting Section and Glossary for future reference.
- Remind participants that the Operator’s Manual is the most comprehensive reference guide for the 2008T BlueStar Machine.

Distribute a *2008T Hemodialysis Machine Quick Start Guide* (P/N 490161) to each participant

- Share that many of the key demonstrations during this training are detailed in the Quick Start Guide (QSG); it is designed to be a ready job aid, but does not substitute for the Operator’s Manual.
- Suggest that participants keep the QSG handy for easy reference as they become familiar with the 2008T BlueStar Hemodialysis Machine.

Notes:

Review 2008T BlueStar Hemodialysis Machine



Review the physical features of the 2008T BlueStar machine, being sure to point out specific features as detailed below:

Back of the machine:

- Power switch
- Heater switch
- 9-volt battery
- USB
- RS232 port
- Incoming water line
- Drain line
- Fresh dialysate line
- Spent dialysate line
- DIASAFE® Plus filter/cover

Right side of the machine:

- Adjustable IV pole height
- Priming bucket with line guides
- Shunt interlock
- Cap holder
- PatientCard Reader
- Incoming dialysate line (BLUE): Sample port and flow meter
- Outgoing dialysate line (RED)
- Brake pedal

Left side of machine:

- Blood pressure cuff basket

Front of machine:

- Display screen
- Status indicator light: **Green** - normal; **Yellow** - warning; **Red** - alarm
- Control Panel Keys: Power; Mute; Reset; Stat Deflate; New Treatment Prime, UF On/Off







Note: Be sure to reinforce the proper sequence for pressing Control Panel keys, refer to the Overview and Setting Parameters sections of the Operator's Manual.

- Touchpad: Confirm; Escape
- Keyboard: Key placement; 10-key panel; Toggle switch for CDX Module
- Modules: Arterial Chamber Holder; Line Guides; Blood Pump; Speed Arrows; Blood Flow Rate (BFR) Display; Pump Segment Settings; Aterial TP Port; Heparin Pump; Air Foam Detector: Air Sensor Heads, Venous Line Clamp, Optical Detector, Line Guides, Venous TP Port.
- Acid connector: O-ring
- Bicarb connector: O-ring
- Ultrafiltration: Fluid sample port



- Encourage participants to take notes on the diagrams contained in the *Participant Workbook*.

Daily Preparation for Treatment

	<p>Refer participants to the detailed instructions found in the Quick Start Guide (QSG). The QSG will be a good resource for participants following the training, so encourage participants to watch the demonstration and use their Participant Workbook to take any additional notes.</p>
	<p>Using the Step-by-Step Instructions Found in the Quick Start Guide:</p>
	<ul style="list-style-type: none"> • Demonstrate how to set up the concentrate • Demonstrate how to set up the bloodline <ul style="list-style-type: none"> o Proper loading and unloading of the blood pump segment o Verify blood pump segment millimeter size • Demonstrate how to prime the extracorporeal circuit <ul style="list-style-type: none"> o Manual or Auto Prime o Remind participants to ensure all air is out of the blood pump segment prior to connecting the arterial line to the dialyzer
	<ul style="list-style-type: none"> • Demonstrate how to test the machine <ul style="list-style-type: none"> o Manual or Auto Start o Remind participants of the criteria for testing: TPs off the machine, stable conductivity, dialysate connector on machine, and no alarms
	<ul style="list-style-type: none"> • Demonstrate the recirculation procedure • Return demonstration is recommended for each of the above steps, depending on the level of experience of the participants. Return demonstrations may be done in either the training area or the treatment setting, at the discretion of the unit.
	<p>Point out topics in the Quick Start Guide to encourage participants to become familiar with using the it as a reference.</p>

Notes:

Review Screens



Review each screen, providing appropriate detail regarding the selections and actions that are available for each screen.

Select Program Screen

- Dialysis
- Rinse
- Acid Clean
- Heat Disinfect
- Chemical Rinse
- Chemical Dwell
- Last Disinfection Cycle and Time
- Software Version
- SLED
- Disinfect Log
- Serial Number
- Apps Installed
- Hour Meter
- PM Due
- DIASAFE Plus Maintenance Due

Dialysate Screen

- Changing Concentrate
- Na+ Level
- Bicarb Level
- Acid/Bicarb Alert
- Alarm Width
- Alarm Position
- Lock/Unlock
- Actual vs. Theoretical Conductivity (TCD)
- Granuflo® Option
- Citrasate Option

Home Screen

- Arterial Pressure Display
- Venous Pressure Display
- TMP Display
- Treatment Clock
- BFR
- Ultrafiltration (UF) Goal
- UF Time
- UF Rate
- UF Removed
- UF Profiles: 1-8
- Dialysate Flow: Auto Flow (All 3) Idle Mode, Sequential
- Temperature
- Conductivity
- Remaining Time of Dialysis (RTD)
- Alarm Volume Control
- Sodium Variation System (SVS) Profiles: Step, Linear, Exponential, None, Starting Na+, Ending Na+, SVS Time



Trends Screen

- Arterial, Venous, and TMP Readings
- Adequacy Monitoring
- SVS/Ultrafiltration Graphics
- Blood Pressure Graphic
- Current vs. Previous Patient Data

Test/Options Screen

- Test: Both Test
- Alarm Test (including 9V battery)
- Pressure Test
- DIASAFE Plus Test
- Time Stamp
- Names of Tests Blue vs. Red "x" as Test Completes
- Independent Conductivity
- Options: Patient ID
- Recirculate
- High Flux
- Single Needle
- Recirculate
- Arterial Alarm Limits
- Venous Alarm Limits (Locked vs. Changeable Limits)
- Verify Crit-Line
- Low Volume
- Auto Prime
- Assisted Reinfusion

Review Screens (continued)

	Heparin Screen	<ul style="list-style-type: none"> • Syringe Type • Hourly Rate • Infusion Time 	<ul style="list-style-type: none"> • Total Infusion • Bolus Amount • Infuse Bolus 	<ul style="list-style-type: none"> • Heparin Prime • Load Syringe • Heparin Dwell (if Selected)
	Kt/V/AF Screen	<ul style="list-style-type: none"> • Online Clearance (OLC) Volume OLC Data Sub-screen: • Number of Tests • Blood Volume Processed (BVP) • Target Kt/V • Delivered Kt/V • Projected Kt/V • Graphic Display: Red vs. Blue Line 	<ul style="list-style-type: none"> • Volume • Actual Time • Plasma Na+ • Yes/No • BVP • Target Kt/V 	<ul style="list-style-type: none"> • Delivered (sp Kt/V) • Delivered (E Kt/V) • Delivered Kt/V • Access Flow
	Blood Volume Monitor (BVM)/Blood Temperature Monitor (BTM) Screen	<ul style="list-style-type: none"> • BVM: Relative Blood Volume (RBV), Arrow Direction Relative Hemoglobin (Hgb)/Hematocrit (Hct) Alert Level, Graphic Display Biofeedback 	<ul style="list-style-type: none"> • BTM: Thermal Energy vs. Temperature Control Recirculation Value Graphic Display 	
	Crit-Line®	<ul style="list-style-type: none"> • Elapsed Time • Initial Hct • Estimated Hgb • Profile; Letter • Current Hct • Current Hgb 	<ul style="list-style-type: none"> • Blood Volume % • Oxygen Saturation • Oxygen Minimum Saturation • Treatment Clock; BV Alert Level • O₂ Alert Level • Print/Markers Button 	<ul style="list-style-type: none"> • BV Graph • O₂/BP Graph • Home Screen • Treatment Information
	Blood Pressure Screen	<ul style="list-style-type: none"> • Systolic, Diastolic, Pulse Upper and Lower Limits • Inflation Pressure Setting 	<ul style="list-style-type: none"> • Clock/Interval Setting • Machine Time Setting/Change 	<ul style="list-style-type: none"> • Graphic Display During Treatment • Treatment Start Time
	<p>Time permitting, ask participants to demonstrate entering treatment parameters and confirming settings after the training is complete. Otherwise, this can be done in the treatment area.</p>			
	<p>Many of these screens may be discussed during the Daily Preparation for Treatment, depending on the trainer. Ensure that all screens have been covered during the In-service.</p>			

Review Alarms: Troubleshooting Activity









Since it is often difficult to simulate alarms and troubleshooting during a demonstration, the following activity has been designed to engage participants in a discussion of alarms. It also helps to reinforce the use of the Troubleshooting Guide and the Operator's Manual to troubleshoot common issues that may occur during treatment.

- Review the reset button; Demonstrate how to reset or hold to completely clear an alarm.
- Introduce the 2008T Troubleshooting Guide and briefly review the topics.
- Run Troubleshooting Activity:
 - o Distribute the guide to the participants. Each participant may have several guides, depending on the size of the group.
 - o Call out an alarm and ask the participant with the guide that matches that alarm to read the information on the guide for the group.
 - o Check to ensure that the answer was clear before moving on to the next alarm.
 - o Address 4-5 key alarms in the same manner, and cover at least one:
 - Blood-related alarm
 - Dialysate-related alarm
 - Warning message
- To close the activity, reinforce that the *Troubleshooting Guide* is designed to be a quick reference, but the *Operator's Manual* is the most comprehensive resource for troubleshooting.
- Collect the *Troubleshooting Guide* and leave it with the unit.
- Continue working with alarms and *Troubleshooting Guide* in treatment area as alarms are encountered (i.e., during return demonstrations).



- This is an excellent opportunity to engage participants; try to select alarms to troubleshoot such that several participants have an opportunity to share with the group.
- Leave the *Troubleshooting Guide* with the unit. If a facility would like additional sets, they may be ordered through the parts department, or the original guide may be copied and laminated by the unit.

Disinfection and Maintenance

    	<ul style="list-style-type: none"> • Using the step-by-step instructions found in the Quick Start Guide, review the procedure for Acid Clean and Heat Disinfect, Auto Heat, if applicable. <ul style="list-style-type: none"> ○ Acid clean is only recommended when closing down the machine for the day. It does not need to be done when bleaching in between patients. ○ Acid clean does not need to occur if bicarbonate was not introduced into the machine. Performing Heat Disinfect and/or Chemical Rinse on spare machines or machines not used, is all that is recommended. • Using the step-by-step instructions found in the Quick Start Guide, review the procedure for Acid Clean and Chemical Rinse (bleach). <ul style="list-style-type: none"> ○ Introduce DIASAFE Plus Filter testing. • DIASAFE Plus filter replacement <ul style="list-style-type: none"> ○ Determine appropriate unit personnel with responsibility for changing filter, DIASAFE Plus Maintenance Reminder, maintenance interval, and DIASAFE Plus reset. ○ Demonstrate removing and replacing the filter per the <i>Operator's Manual</i>. <ul style="list-style-type: none"> - The filter should be disposed of per the unit's policy. - Heat Disinfect should be run after replacing the filter to disinfect the machine. ○ Demonstrate testing the DIASAFE Plus filter. <ul style="list-style-type: none"> - Note: The DIASAFE Plus filter test shows a red "X" when setting up the machine, the operator should perform a test before moving forward to ensure that the filter is functioning. - Remind participants that machines should never be used for patients if there are red "X"'s in the test and option screen. • Return demonstration is recommended either at the end of the In-Service, or in the treatment area at the end of the day.
	<ul style="list-style-type: none"> • Because responsibility for cleaning and maintenance varies by facility, be sure to check with each unit to determine appropriate personnel to train, as well as policies and procedures.

2008T BlueStar Hemodialysis Machine In-Service Training Checklist



Facility: _____ Date of In-Service Training: _____

Address: _____ Contact: _____

City: _____ State: _____ Zip: _____ Phone: _____

Presented by: _____, RN

Training Topics

Topics may be modified to accommodate the specific policies and procedures of the unit being trained. Check marks indicate topics reviewed during the In-Service training; no representation is made regarding the participant competence or expertise.

Training Topics	✓		✓
1. 2008T BlueStar Features			
• Familiarize Physical Features		• Low Power Mode	
• Review of Screen Attributes		• Idle Mode Option of 100 mL/min. Dialysate Rate	
• Review of Alarms/Troubleshooting		• New 150 mL/min. Dialysate Rate	
• Disinfection & Maintenance		• Isolated UF	
• Auto Start/Testing		• Alarm Volume Control, includes: Difference Between Low and High Priority Sounds, Acute vs. Standard (Chronic), How to Adjust Volume, Mute Symbol, and Volume Control Limits.	
• Sustained Low Efficiency Dialysis (SLED)			
• Low Volume Mode		2. Service Mode Options	
• Disinfect Log		• Default Dialysate Flow 500/600/700/800 mL/min.	
• Auto Prime		• Plasma Sodium — Hide/Show	
• Assisted Reinfusion		• Sodium Variation System (SVS) On — Yes/No	
• Independent Conductivity/pH Testing		• Idle Mode Dialysate Flow Rate — 100/300 mL/min.	
• PatientCard System		• Conductivity Alarm Position — Lock/Unlock	
– Entering Patient Information		• New Treatment Prescription Warning	
– Default Parameters		• Default Rx Screen	
– Dialysate/Heparin/BP/UF Settings		• Low Volume Max UF Rate — 500 – 1,000 mL	
– Prescription		• Maximum UF Per Treatment vs. Per Hour	
• Alarm Volume Controls		• Alarm Type — Acute/Chronic	
• Select Program Screen Changes		• Maintenance Scheduler Options	
– DIASAFE® Plus Reminder			
– PM Reminder			
– Digital Hour Meter			
– Digital Serial Number			
• New Heparin Syringe Choices			
• Conductivity Alarm Position Lock (if clinic is using)			
• New Treatment Prescription Warnings			



**FRESENIUS
MEDICAL CARE**

RENAL TECHNOLOGIES

Fresenius Renal Technologies, a division of Fresenius Medical Care North America
920 Winter Street, Waltham, MA 02451
Customer Service: 800-662-1237 | Technical Support: 800-227-2572

2008TBlueStar.com