Akut-Therapie-Systeme

Ultraflux EMI[®]2

Erhöhte Mittelmolekül Clearance

 multi**Intense**Care

**FRESENIUS
MEDICAL CARE**

Inhalt

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01 Susantitaphong P et al. World Incidence of AKI: A Meta-Analysis. *Clin J Am Soc Nephrol* 2013;8(9):1482-1493. **02** Hoste et al. Epidemiology of acute kidney injury in critically ill patients: the multinational AKI-EPI study. *Intensive Care Med.* 2015;41(8):1411-1423. **03** Siebeck M et al. Super high-flux CVVHD using regional citrate anticoagulation: long-term stability of middle molecule clearance. *Crit Care.* 2015;19(Suppl 1):P301. **04** Morgera S et al. Renal Replacement Therapy With High-Cutoff Hemofilters: Impact of Convection and Diffusion on Cytokine Clearances and Protein Status. *AJKD.* 2004;43(3):444-453. **05** Kidney Disease: Improving Global Outcomes (KDIGO) Acute Kidney Injury Work Group. KDIGO Clinical Practice Guideline for Acute Kidney Injury. *Kidney Int Suppl.* 2012;2:1-138. **06** Wald R et al. The Association Between Renal Replacement Therapy Modality and Long-Term Outcomes Among Critically Ill Adults With Acute Kidney Injury: A Retrospective Cohort Study. *Crit Care.* 2014;42(4):868-877. **07** Vanholder R et al. Rhabdomyolysis. *J Am Soc Nephrol.* 2000;11:1553-1561. **08** Fabbri P et al. Myoglobin removal of small-protein leakage membrane (EMIC2) in patients in the ICU: a case series. *Crit Care.* 2014;18(Suppl 1):P398. **09** Honoré PM et al. Prevention and treatment of sepsis-induced acute kidney injury: an update. *Ann Intensive Care.* 2015;5(51):1-10. **10** Rimmelé T et al. Super high-flux continuous hemodialysis: an efficient compromise for blood purification in sepsis. *Crit Care.* 2012;16(Suppl 1):P377. **11** Chatzizisis YS et al. The syndrome of rhabdomyolysis: complications and treatment. *Eur J Intern Med.* 2008;19(8):568-574. **12** Bellomo R et al. Acute kidney injury in sepsis. *Intensive Care Med.* 2017;43(6):816-828. **13** Sancharawala V. Light-Chain (AL) Amyloidosis: Diagnosis and Treatment. *Clin J Am Soc Nephrol.* 2006;1:1331-1341. **14** Bellomo R et al. Renal replacement therapy in the intensive care unit. *Crit Care and Resuscitation.* 1999;1:13-24. **15** Bellomo R et al. Nomenclature for continuous renal replacement therapies. *AJKD.* 1996;28(5):S1-S7. **16** Honoré, et al. The big bang of hemofiltration: The beginning of a new era in the third millenium for extra-corporeal blood purification. *Int J Artif Organs.* 2006;29(7):649-659. **17** Schmidt JJ et al. New high-cut-off dialyzer allows improved middle molecule clearance without an increase in albumin loss: a clinical crossover comparison in extended dialysis. *Blood Purif.* 2012;34:246-252. **18** Kuzmanovska B et al. Rhabdomyolysis in critically ill surgical patients. *Med Arch.* 2016;70(4):308-310. **19** Guzman N et al. Myoglobin removal using high-volume high-flux hemofiltration in patients with oliguric acute kidney injury. *Blood Purif.* 2013;36:107-111. **20** Bosch et al. Rhabdomyolysis and acute kidney injury. *N Engl J Med.* 2009;361:62-72. **21** De Meijer AR et al. Serum creatine kinase as predictor of clinical course in rhabdomyolysis: a 5-year intensive care survey. *Intensive Care Med.* 2003;29:1121-1125. **22** Petejova N et al. Acute kidney injury due to rhabdomyolysis and renal replacement therapy: a critical review. *Crit Care* 2014;18:224. **23** Sorrentino SA et al. High permeability dialysis membrane allows effective removal of myoglobin in acute kidney injury resulting from rhabdomyolysis. *Crit Care Med.* 2011;39(1):184-186. **24** Leppänen I et al. Dialyzer clearance of myoglobin with middle molecule filter and low blood flow CVVHD in patients with rhabdomyolysis-associated acute kidney injury. *Crit Care.* 2013;17(Suppl 2):P432. **25** Cantaluppi V et al. Protective effect of high molecular flow membranes on inflammatory and myoglobin-associated acute kidney injury in patients with rhabdomyolysis after major cardiac surgery. *NDT.* 2015;30(Suppl 3):iii133-iii134. **26** Singer M et al. The Third International Consensus Definitions for Sepsis and Septic Shock (Sepsis-3). *JAMA.* 2016;315(8):801-810. **27** Kellum JA et al. Understanding the inflammatory cytokine response in pneumonia and sepsis. *Arch Intern Med.* 2007;167(15):1655-1663. **28** Ronco C et al. Interpreting the Mechanisms of Continuous Renal Replacement Therapy in Sepsis: The Peak Concentration Hypothesis. *Artificial Organs.* 2003;27(9):792-801. **29** Engel et al. Epidemiology of sepsis in Germany: results from a national prospective multicenter study. *Intensive Care Med.* 2007;33(4):606-618. **30** Gordon AC et al. Effect of Early Vasopressin vs Norepinephrine on Kidney Failure in Patients With Septic Shock The VANISH Randomized Clinical Trial. *JAMA.* 2016;316(5):509-518. **31** Rhodes A et al. Surviving Sepsis Campaign: International Guidelines for Management of Sepsis and Septic Shock: 2016. *Intensive Care Med.* 2017;43(3):304-377. **32** Eichhorn T et al. Clearance of selected plasma cytokines with continuous venovenous hemodialysis using Ultraflux EMIC2 versus Ultraflux AV1000S. *Blood Purif.* 2017;44:260-266. **33** Kalb R et al. Regional citrate anticoagulation for high volume continuous venovenous hemodialysis in surgical patients with high bleeding risk. *Ther Apher Dial.* 2013;17(2):202-212. **34** Anderson KC et al. NCCN Guidelines Insights. Multiple Myeloma. Version 3.2016. *J Natl Compr Canc Netw.* 2016;14(4):389-400. **35** Singhal S et al. Multiple myeloma. *Clin J Am Soc Nephrol.* 2006;1(6):1322-1330. **36** Hanf W et al. Prolonged hemodialysis for acute kidney injury in myeloma patients. *Clin Nephrol.* 2009;74(4):319-322. **37** Falk RH. Diagnosis and management of the cardiac amyloidoses. *Circulation.* 2005;116(20):2047-2060. **38** Knudsen LM et al. Renal failure in multiple myeloma: reversibility and impact on the prognosis. *European J Haem.* 2000;65(3):175-181. **39** Bladé J et al. Renal failure in multiple myeloma: presenting features and predictors of outcome in 94 patients from a single institution. *Arch Intern Med.* 1998;158(17):1889-1893. **40** Zannetti BA et al. Bortezomib-based therapy combined with high cut-off hemodialysis is highly effective in newly diagnosed multiple myeloma patients with severe renal impairment. *Am J Hematol.* 2015;90(7):647-652. **41** Hutchison CA et al. Efficient Removal of Immunoglobulin Free Light Chains by Hemodialysis for Multiple Myeloma: In Vitro and In Vivo Studies. *J Am Soc Nephrol.* 2007;18(3):886-895. **42** Morgera S et al. A Safe Citrate Anticoagulation Protocol with Variable Treatment Efficacy and Excellent Control of the Acid-Base Status; *Crit Care Med.* 2009;37(6):2018-2024. **43** Morgera S et al. Metabolic Complications during Regional Citrate Anticoagulation in Continuous Venovenous Hemodialysis: Single-Center Experience. *Nephron Clin Pract.* 2004;97:c131-c136. **44** Link A et al. Total-to-ionized calcium ratio predicts mortality in continuous renal replacement therapy with citrate anticoagulation in critically ill patients. *Crit Care.* 2012;16:R97. **45** Schultheiss C et al. Continuous venovenous hemodialysis with regional citrate anticoagulation in patients with liver failure: a prospective observational study. *Crit Care.* 2012;16:R162. **46** Raimundo M et al. Maintaining normal levels of ionised calcium during citrate-based renal replacement therapy is associated with stable parathyroid hormone levels. *Nephron Clin Pract.* 2013;124:124-131. **47** Khadzhynov D et al. Incidence and outcome of metabolic disarrangements consistent with citrate accumulation in critically ill patients undergoing continuous venovenous hemodialysis with regional citrate anticoagulation. *J Crit Care.* 2014;29(2):265-271. **48** Tovey L et al. Beyond the randomised clinical trial: citrate for continuous renal replacement therapy in clinical practice. *Nephron Clin Pract.* 2013;124:119-123. **49** Oudemans-van Straaten HM et al. Citrate Anticoagulation for Continuous Renal Replacement Therapy in the Critically Ill. *Blood Purif.* 2010;29:191-196. **50** Wu MY et al. Regional Citrate Versus Heparin Anticoagulation for Continuous Renal Replacement Therapy: A Meta-Analysis of Randomized Controlled Trials. *Am J Kidney Dis.* 2012;59(6):810-818. **51** Bai M et al. Citrate versus heparin anticoagulation for continuous renal replacement therapy: an updated meta-analysis of RCTs. *Intensive Care Med.* 2015;41:2098-2110. **52** Liao YJ et al. Citrate versus unfractionated heparin for anticoagulation in continuous renal replacement therapy. *Chin Med J.* 2013;126(7):1344-1349. **53** Liu C et al. Regional citrate versus heparin anticoagulation for continuous renal replacement therapy in critically ill patients: a meta-analysis with trial sequential analysis of randomized controlled trials. *Crit Care.* 2016;20:144. **54** Instruction for Use, Ultraflux EMIC2, 08/2017. **55** Klingel R et al. Comparative analysis of procoagulatory activity of haemodialysis, haemofiltration and haemodiafiltration with a polysulfone membrane (APS) and with different modes of enoxaparin anticoagulation. *Nephrol Dial Transplant.* 2004;19(1):164-170. **56** Friedrich JO et al. Hemofiltration compared to hemodialysis for acute kidney injury: systemic review and meta-analysis. *Crit Care.* 2012;16:R146. **57** Ricci Z et al. Solute Removal during continuous renal replacement therapy in critically ill patients: convection versus diffusion. *Crit Care.* 2006;10(2):R67.

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Firmensitz: Fresenius Medical Care Deutschland GmbH · 61346 Bad Homburg v. d. H. · Deutschland
Deutschland: Fresenius Medical Care GmbH · Else-Kröner-Straße 1 · 61352 Bad Homburg v. d. H.
Österreich: Fresenius Medical Care Austria GmbH · Lundenburgergasse 5 · 1210 Wien
Schweiz: Fresenius Medical Care (Schweiz) AG · Aawasserstrasse 2 · 6370 Oberdorf NW