

Fresenius Medical Care North America Service Bulletin


Equipment: All Hemodialysis Machines
Bulletin: 14-FRS-002 Rev A
Subject: Air Separator Adapter Board

1.0 PURPOSE

To inform all Field Service and Technical Support personnel of the part number for the Air Separator Adapter Board Kit.

2.0 SCOPE

2008[®] Series Hemodialysis machines that have been produced or retrofitted with the Air Separator Adapter Board will have a small board installed on the machine's distribution board under the air separator connector (x6).



Replacement or spare Air Separator Adapter Boards are available by ordering the Air Separator Adapter Board Kit (P/N 191011). The kit includes an Air Separator Adapter Board and installation instructions.

Contact Fresenius Medical Care Customer Service Spare Parts at 800-227-2572 to order.

If there are any questions regarding this bulletin, contact Fresenius Medical Care Technical Support at 800-227-2572

© 2014, Fresenius Medical Care North America, All rights reserved.