



## Fresenius Medical Care North America Service Bulletin

Equipment: Dry Acid Dissolution Units  
Bulletin: 24-DRY-001 Rev A  
Subject: Surge Protector Recommendation

### 1.0 PURPOSE

To inform all Field Service and Technical Support personnel of a recommended surge protector to be used with the 99 gallon and 132 gallon Dry Acid Dissolution Units (DADU).

### 2.0 DESCRIPTION

If there is any concern that the facility's electrical supply contains the occasional surge or electrical spike, Fresenius Medical Care recommends the use of a surge protection device like the [Tripp Lite Model ISOBAR6ULTRAHG 6-Outlet Surge Protector](#) or equivalent.



Tripp Lite  
Model ISOBAR6ULTRAHG  
6-Outlet Surge Protector

© 2024 Fresenius Medical Care. All Rights Reserved.

If there are any questions regarding this bulletin, contact Fresenius Medical Care Technical Support at 800-227-2572.