

# 6008 CAREset

## Product Range

**All-in-one disposable with completely pre-connected bloodlines for all treatment modalities**

### Air free design and low extracorporeal volume

- Innovative technologies lead to an air-free<sup>1</sup> and flow-optimized design of the extracorporeal system
- Low extracorporeal volume and reduced contact between blood and foreign material
- Both (less blood-air contact and low extracorporeal volume) may reduce the risk of clotting during therapy<sup>2</sup>

### Economic efficiency and eco-friendliness

- All-in-one CAREset simplifies logistics
- Compact design and lightweight material of CAREset reduces waste



reddot design award  
winner 2018<sup>3</sup>



Product	Product code	Inner diameter pump segment	Filling volume (max.)		Units/box
			AV*	SN**	
6008 CAREset-R	F00006808	7.1 mm	115 mL	127 mL	26
6008 CAREset BVM-R	F00007727	7.1 mm	122 mL	134 mL	26
6008 CAREset L-R	F00007728	7.1 mm	128 mL	140 mL	26
6008 CAREset BVM L-R	F00007729	7.1 mm	135 mL	146 mL	26
6008 CAREset Low Volume BVM-R	F00007579	7.1 mm	83 mL		24
6008 CAREset Low Volume-R	F00007578	7.1 mm	78 mL		24

\* Double-needle treatment

\*\* Single-needle treatment, without stroke volume

#### References:

<sup>1</sup> Valid for Double-needle treatment

<sup>2</sup> Kessler M et al. Anticoagulation in Chronic Hemodialysis. Semin Dial. 2015 Sep-Oct; 28(5) 474-89

<sup>3</sup> reddot design award is an international design competition, dating back to 1955. It is one of the most sought-after seals of quality for good design.

