




**reddot** design award  
winner 2018\*



- **Post-dilution port**  
Easy to switch between treatment modalities
- **Integrated infusion port**  
Easy delivery of HighVolumeHDF®
- **Single-Needle chamber**  
Maximum flexibility in the delivery of therapy
- **Pre-dilution port**  
Easy to switch between treatment modalities
- **Venous chamber**  
Innovative technologies for an air-free design
- **Air-free pressure measurement (pre-filter)**  
Reduced risk of blood coagulation
- **Air-free pressure measurement (arterial)**  
Reduced risk of blood coagulation
- **Air-free pressure measurement (venous)**  
Reduced risk of blood coagulation

*\*Information: The Red Dot Design Award is an international design competition, dating back to 1955. It is one of the most sought-after seals of quality for good design. The 6008 CAREset won this award in the "Medical devices and technology" category. For more details, please visit: <https://www.red-dot.org>*