






Accessories

Reliability established over many years
and millions of treatments

Type	Target system	Sterilization	Description	Depiction	Units per box	Article code
Twister	universal	ETO	Enables reversal of arterial and venous needle at patient access to measure fistula flow in combination with access flow measurement software feature		50	F00000193
Recirculation connector	universal	Beta	Adapter to connect with two male Luer-Lock connectors		100	F00007581
Spike	universal	Beta	Spike to connect to an infusion bag		100	F00007582
Vented Spike	universal	ETO	Spike to connect to an infusion bag/bottle		12	F00007583
Drainage Bag 2lt	universal	ETO	Collection bag for priming fluid with female Luer-Lock and clamp		100	5015091
Drainage connector	5008	Beta	Rinse connector for ONLINE priming with 5008 Therapy System		12	7030031
Rinse connector	6008	Beta	Rinse connector for ONLINE priming with 6008 CAREsystem		12	F00009519
SN-Adapter Luer-Lock	universal	ETO	Y-Adapter for single needle treatment with attached tubes and clamps		100	5027851
Frekaflex infusion set	universal	ETO	Infusion system with roller clamp and drip chamber		150	2889011
Dialysis accessories	universal	ETO	Infusion set with drip chamber, roller clamp, spike, enclosed recirculation connector and additional spike		130	5017151
Transducer Protector	universal	ETO	Transducer protector with female Luer-Lock male Luer-Lock connector		100	F00007580
Pressure line for MTS	universal	ETO	Transducer protector with 30 cm tube		100	5014631
Dialyser Exchange-Set	universal	ETO	Tubing system for exchange of dialyser during treatment, short-circuit AV-set		75	5018001
Fluid Sampling Set	5008	Beta	Fluid sampling set for quality testing with 5008 Therapy System		15	F00006222
6008 Fluid Sampling Set	6008	Beta	Fluid sampling set for quality testing with 6008 CAREsystem		15	F00010979
Adapter Male Luer-Lock/ Male Luer-Lock	universal	ETO	Adapter to connect with two female Luer-Lock connectors		100	5014771
Prolongation LLF-LLM 75 cm	universal	Beta	Prolongation tube with female Luer-Lock and male Luer-Lock connector		100	7030011