



Cartridge-172 Identification Tool



1	Post-filter Port (blue clamp)	6	SecureClip	11	Arterial Patient Line (red clamp)	16	Dialysate Inlet Line (green clamp)
2	Venous Needleless Access Port	7	Female-female Connector	12	Saline "T" (white clamp)	17	Effluent Line
3	Heparin Line	8	Saline Line (white clamp)	13	Waste Line (yellow clamp)	18	Cartridge Label
4	Dialysate Outlet	9	Priming Spike	14	Saline Needleless Access Port		
5	Venous Patient Line (blue clamp)	10	Connector Tab	15	Access Pressure Pod		

Cartridge-172 Identification Tool





Cartridge-172 Identification Tool

Priming Spike	Connector Tab	Arterial Patient Line (red clamp)	Saline “T” (white clamp)
Female – female connector	Access Pressure Pod	Arterial Air Detector	Effluent Line
Post-dialyzer Port (blue clamp)	Cartridge Information Label	Dialysate Fluid Outlet Line	Venous Air Detector
Dialysate Fluid Air Detector	Dialysate Fluid Inlet Line (green clamp)	Saline Line (white clamp)	Venous Patient Line (blue clamp)
Connector Tab	Waste Line (yellow clamp)	Venous Needleless Access Port	Saline Needleless Access Port
Heparin Line	SecureClip		

(Labels formatted for Avery Label 1” x 2.63”)