



NxRx®

Online Prescription System



Introduction to NxRx	Transactions	Nx2me®
Get Help	Reports	Contact Management
Dashboard	Documents	My Account
Patients	Move Equipment	

Introduction to NxRx

NxRx Key Features

App Update Notification

Overview

Signature Expiration

Logging In

Layout Overview

Terms and Conditions



Overview

NxRx is a comprehensive online platform that supports prescription and operational management for home hemodialysis (HHD) therapy.

Designed for use with NxStage® systems, NxRx enables healthcare providers to create, manage, and adjust prescriptions while streamlining supply fulfillment for home patients. Beyond prescription management, the platform offers tools for generating reports, managing documents, and tracking center inventory—all in one centralized system

Sample screenshots are for illustrative purposes only and do not reflect real patient or medically relevant data.



NxRx Key Features

Patient and Prescription Overview:

Displays detailed patient information alongside their associated prescriptions

Form Automation:

Automatically pre-fills forms using existing patient and prescription data

Web Accessibility:

Web-based platform accessible from any PC with an internet connection

Comprehensive Reporting Modules:

- Equipment management records
- Patient compliance metrics
- Patient therapy configuration details
- Dialysate bag usage statistics
- Quality assurance release certification



Logging In

The login form is titled "Welcome to NxRx". It contains four input fields, each highlighted with a red border: "Account ID", "Username", "Password", and a "Log In" button. The "Password" field has a "Show" button to its right. Below the "Log In" button is a link "Forgot your password?". At the bottom of the form is a link "Need Help Logging In?".

- Go to <https://my.nxstage.com> in your web browser
- Enter the following:
 - Your NxStage® assigned facility **Account ID**
 - Your email address as your **Username**
 - Your NxStage assigned **password**
- Click **Log In**

For access or assistance, contact Customer Service at **866-697-8243** (Option **2**) or customerservice@fmc-na.com



Terms and Conditions

Terms and Conditions

Welcome to NxStage Medical's electronic prescription section. From here you'll be able to enter prescription information for your patients. This section of our site is only available to dialysis centers that have an authorized physician signature on file with us.

By selecting the agree button below you acknowledge that an authorized physician signature is on file and agree to the terms and conditions of use below.

This system is restricted to authorized users for legitimate business purposes only. By accessing this site you agree to use the NxStage website only for lawful purposes. Your password is intended to protect your organizations' privacy and should be kept confidential.

Unauthorized attempts to access, use or modify wvm.nxstage.com content is a violation of federal and state laws.

If you have any questions, please contact our Customer Service department at 1-866-NXSTAGE (697-8243).


☒ I Agree
 [Cancel](#)



- Review the Terms and Conditions
- Check the **I Agree** box
- Click **Accept**




App Update Notification



- Review any changes or updates to the application
- Click **Continue**



Signature Expiration Notification



- Review any physician signatures that are expired or nearing expiration
- Click **Okay**



Layout Overview

St. Elmo (MA)
438 South Production St.
Kansas City, Massachusetts 12345
(518) 453-1669 x4545

A

Log Out

B

Get Help | Dashboard | Patients | Welcome Calls | Transactions | Reports | Documents | Move Equipment | Nx2me | Contact Management | My Account

C

New Patient | New ETD Patient | Show/Hide Dropped Patients | Bag Approval 0

D

Search

NxStage ID	Patient Name	DOB	Status	Start Date	Valid Thru	Action
39638.0	Demo, Frank	1/1/1950	Active	5/1/2017	Until Notified	View Current Rx Revise Change Status Travel
99002.0	DemoPatient, Amy	4/1/1948	Active	4/1/2014	Until Notified	View Current Rx Revise Change Status Travel
99001.0	DemoPatient, Charles	4/1/1948	Active	4/1/2014	Until Notified	View Current Rx Revise Change Status Travel
125000.0	DemoPatient, Christopher		Active			View Current Rx Revise Renew Change Status Travel

E

Feedback

G

To access the NxStage Dosing Calculator, click [here](#)

Invent. Improve. Inspire.

- A. Log Out
- B. Main menu
- C. Action buttons
- D. Search bar
- E. Nx2me App Feedback
- F. Prescription Quick links
- G. Dosing Calculator link



Get Help

Get Help



Get Help

Subject: Help me with NxRx Online Orders

From: Peter North

Phone:

Description:

You have used 0 of the possible 255 characters that are allowed.

 Submit  Cancel

- Provide your phone number and a brief description of your needs
- Click **Submit**
- Receive a response within 24 hours



Dashboard

Dashboard



- Get Help
- Dashboard
- Patients
- Welcome Calls
- Transactions
- Reports
- Documents
- Move Equipment
- Nx2me
- Contact Management
- My Account



Dashboard Actions

Bag Approvals 0	Expiring Physician Signatures 0	Expired Prescriptions 5	Pending RMAs 0	On Hold > 90 Days 0
Unemployed Equipment 0	Unsubmitted Transactions 57	Upcoming Service Swaps 0	Upcoming Travel Orders 1	Signature Required 39

- Review the dashboard for important actions, including expired prescriptions and signatures required
- Click a dashboard tile to open the page containing the displayed data



Patients

Patient Overview

Adding a New Patient

Prescription Quick Links

Bag Approvals

Finding unlisted patients



A

New Patient

B

Show/Hide Dropped Patients

C

Bag Approval 0

NxStage ID	Patient Name	DOB	Status	Start Date	Valid Thru	Action
Search D						
11802.3	Bush, Petunia	10/30/1980	Active	12/9/2023	12/31/2025	View Current Rx Revise Change Status Travel
99999.6	Bush, Petunia	10/30/2012			1/25/2023	View Current Rx Add
99002.0	DemoPatient, Amy	4/1/1950	Active	4/1/2014	Until Notified	View Current Rx Revise Change Status Travel
99001.0	DemoPatient, Charles	4/1/1948	Active	4/1/2014	Until Notified	View Current Rx Revise Change Status Travel
125000.0	DemoPatient, Christopher		Active			View Current Rx Revise Renew Change Status Travel

- A. Add a prescription for a new patient (Buttons shown based on center configuration)
- B. Toggle **Dropped Patients** visibility
- C. Display **Bag** requests and approvals
- D. **Prescription Quick Links**



Get Help	Dashboard	Patients	Welcome Calls	Transactions	Reports	Documents	Move Equipment	Nx2me	Contact Management	My Account
----------	-----------	----------	---------------	--------------	---------	-----------	----------------	-------	--------------------	------------



[View Current RX](#)

View a patient's current prescription

[Revise](#)

Revise a patient's prescription (add *Transition to Home* date when known)

[Renew](#)

Renew a patient's expired prescription

[Change Status](#)

Notify NxStage of patient hospitalizations, discontinuations, and credit requests

[Travel](#)

Submit a travel prescription form



Enter Prescription

Patient: **Robert Smith**

Patient P.O. & Delivery Support Physician Treatment Equipment Training Supplies Notes

Patient Information

Patient ID: **NEW**

First Name: * Robert Middle: Last Name: * Smith

Date of Birth: * 7/4/1976

Address Validated: **No**

Note: If the address entered is not validated, there is a possibility that the delivery may not reach the correct location.

Address: * 2 monument sq (eg: 3511 Oak Road)

- Enter patient information
- Select **Next**
- Choose **Validate** or **By-Pass Validation**

Address Not Validated ×

The Address entered is not yet validated. Do you want to validate?

Validate **By-Pass Validation**

- **Save** may be selected at any step





Enter Prescription

Patient: Robert Smith

Patient P.O. & Delivery Support Physician Treatment Equipment Training Supplies Notes

Purchase Order and Dates

Purchase Order #:

Center Training Start Date: 

Order Valid Until Otherwise Notified: ☒ Yes ☐ No

Date Shipment Needed: 

Transition To Home: 

- Complete Purchase Order and date section when ordering equipment (Purchase order number required by some centers)
- Select [Welcome Home](#) option
- Enter delivery information and alternate shipping address as needed
- Select **Next**





Get Help	Dashboard	Patients	Welcome Calls	Transactions	Reports	Documents	Move Equipment	Nx2me	Contact Management	My Account
----------	-----------	----------	---------------	--------------	---------	-----------	----------------	-------	--------------------	------------



Welcome Home Call

Call or Video: ☒ Video ☐ Call

Email Address:

This email does NOT need to be the patient's email, however, the link sent to this email will be patient specific.

Welcome Home Call

Call or Video: ☐ Video ☒ Call

Please call NxStage Customer Service during regular business hours 2 weeks prior to your patient transitioning home for an on the spot Welcome Home Call.

Select a *Welcome Home* option:

- **Video** - Enter email address (does not need to be the patient's email, but the generated link will be patient specific)
- **Call** - Call NxStage Customer Service 2 weeks prior to patients transitioning home for a *Welcome Home* call to be initiated
- **Save** may be selected at any step





Enter Prescription

Patient: Robert Smith

Patient P.O. & Delivery Support Physician Treatment Equipment Training Supplies Notes

Therapy Coordinator *

	Last Name	First Name	Position	
☐	Patel	Aditya	Administrator	
☒	Clawson	Alexander	Nurse	(518) 453-1
☐	tejeda	arbelnellys	Home Training Nurse	(518) 453-
☐	Ella	Barbar	Home Training Manager	

- Check the box for the Therapy Coordinator that is entering the Prescription information or is assigned to the patient
- Select Next





Enter Prescription

Patient: **Robert Smith**

[Patient](#) [P.O. & Delivery](#) [Support](#) [Physician](#) [Treatment](#) [Equipment](#) [Training Supplies](#) [Notes](#)

Physician *

	Name	Position	Phone
☐	Roberts, Matthew	Physician	
☐	Breen DO NOT USE, Cormac INACTIVE	Physician	02071885665
☒	Roberts, Matthew	Physician	
☐	Roberts, Matthew	Physician	

- Check the box for the physician associated with the patient. If a signature is required, the system will indicate whether one is on file and display its expiration date
- Select **Next**





Enter Prescription

Patient: **Robert Smith**

[Patient](#) [P.O. & Delivery](#) [Support](#) [Physician](#) [Treatment](#) [Equipment](#) [Training Supplies](#) [Notes](#)

Monthly Treatment Package (Includes up to 26 treatments)

Therapy Type: * ☒ Pureflow SL (Concentrate) ☐ Express (Premixed Bags)

Treatment Schedule: *

Fluid Volume: *

Hanging Bag Volume (Non-PFSL): * ?

Dialysate flow rates greater than 12 Liters (L) per hour (hr): * ☐ Yes ☒ No

- Enter patient treatment details
- Order accessory supplies as needed
- Select **Next**





Enter Prescription

Patient: **Robert Smith**

[Patient](#) [P.O. & Delivery](#) [Support](#) [Physician](#) [Treatment](#) [Equipment](#) [Training Supplies](#) [Notes](#)

NxStage System One Equipment

Source Water Quality Total Dissolved Solids: ppm Hardness: grain

Cycler * ☒ New ☐ Redeployed ☐ N/A (Serial #)

PureFlow™ SL Control Unit * ☒ New ☐ Redeployed ☐ N/A (Serial # D)

- Enter **Source Water Quality** (required for PureFlow patients)
- Order equipment (provide serial numbers for redeployed equipment)
- Select **Next**





Enter Prescription

Patient: **Robert Smith**

[Patient](#) [P.O. & Delivery](#) [Support](#) [Physician](#) [Treatment](#) [Equipment](#) [Training Supplies](#) [Notes](#)

Training Supplies

The Training Supplies tab in NxRx has been created to assist you with the ordering of training supplies. The Accessories and User's Guide section has been pre-populated with all recommended supplies selected. Should you not need one of these items please unselect them.

Accessories and User Guides

Part Number	Description	UoM	Requested
FW-300-1	Express Fluid Warmer Accessories	ea	☒

- Review pre-selected training supplies and edit selection as needed
- Select **Next**





Enter Prescription

Patient: **Robert Smith**

Patient P.O. & Delivery Support Physician Treatment Equipment Training Supplies **Notes**

Enter your notes and comments here:

You have used 0 of the possible 255 characters that are allowed.

← Previous **Save** **Submit**

By his/her signature on file, you (the Prescriber) certify that before NxStage is directed to ship products to the patient's home, the qualified, as required by the System One User Guide and that all disclosures of patient information to NxStage are made in accordance with applicable federal or state law, the use by and disclosure to NxStage of all the information provided in this Additional Patient Information is not required for order fulfillment. By signing below, you (the Prescriber) also certify that you are not providing information to NxStage for the purpose of marketing products, services, and programs, and for market research and analysis.

★ Prescriber Signature

- Enter information advising NxStage of any special processing requirements or additional information
- Select **Save**
- Select **Submit**
- Select **Yes** to confirm
- Your browser downloads a PDF copy
- Select **Okay** to return to the app



Patient Overview

New Patient

New ETD Patient

Show/Hide Dropped Patients

Bag Approval 1

Open Bag Requests

Last Date Bank Updated: 1/31/2016
(Bank balance is updated within the first 10 days of each month)
Bank Amount: 5

Important – Bag authorizations are needed by the end of the next business day. Bags will not be shipped if

Approve

Patient	Number of Cases On Hand	Number of Cases Re
DemoPatient, Elizabeth	0	
		Total Cases Re

Click patients name to view patient's recent bag approval history and current reliability bag activity.

- When there are open bag requests, the **Bag Approval** button turns orange and displays the number of requests
- Click **Bag Approval**
- Select a patient and click **Approve**





Bag Approval

Patient First Name:

Elizabeth

Patient Last Name:

DemoPatient

Number of Cases on Hand:

0

Number of Cases Requested:

15

Number of Reliability Cases That Will Be Shipped Due To A Previous Reliability Event:

0

Notes/Comments From Patients:

Patient's Reason:

Travel

Bags Approved:

Note: If you choose to deny all bags, enter a quantity of zero and press submit.

Submit

Patient Inventory History (Past 90 Days)

Patient	Month Year	Patient's Reason	Number of Cases Shipped
---------	------------	------------------	-------------------------

- Enter the number bags to approve (to deny all bags, enter a quantity of zero)
- Click **Submit**



Standard New Patients

- Under the **Patients** tab, a patient may not be visible until their prescription has been processed which can take up to 24 hours.
- In Salesforce, search for a patient's name or prescription to view the current status.
- **Note:** Prescriptions won't be visible until the **Therapy Start Date**

Transfer Patients

- Check the pending prescription in Matt Robert's Prescription file in Salesforce. The doctors' signature is needed to process the transfer.
- For *transfer* questions, contact Matthew.Roberts@freseniusmedicalcare.com
- **Note:** Transfers cannot be processed until the **Therapy Start Date**



Transactions

Transactions



Review Transactions - Already Submitted

[Go To Not Yet Submitted](#)

[View Transaction](#)

1 to 20 of 1007

Submitted Date/Time	Created By	NxStage ID	Patient Name	Transaction Type	Confirmation Number
5/21/2025, 11:58 AM		11802.3	Petunia Bush	Revise	RX-05212025-492747
4/30/2025, 11:01 AM	Bart Simpson	11802.3	Petunia Bush	Revise	RX-04302025-554801
4/24/2025, 11:30 AM	Bart Simpson	New	Sarah Smith	New Prescription	RX-04242025-554134
4/24/2025, 10:48 AM	Nurse1 ReadOnly	11802.3	Petunia Bush	Change Status: Credit	st-04-24-2025-723638
4/24/2025, 10:47 AM	Nurse1 ReadOnly	11802.3	Petunia Bush	Change Status: Resu...	st-04-24-2025-723637
4/24/2025, 10:46 AM	Nurse1 ReadOnly	11802.3	Petunia Bush	Change Status: Hold ...	st-04-24-2025-723636
4/24/2025, 10:39 AM	Nurse1 ReadOnly	New	April May	New Prescription	RX-04242025-554102

- Toggle between **Already submitted** and **Not yet submitted**
- Select an entry from the list and click **View Transaction**
- **Not yet submitted** transactions also have an option to **Delete**
- **Already submitted** transaction PDFs are available by clicking **Get PDF**



Reports

Reports



Get Help	Dashboard	Patients	Welcome Calls	Transactions	Reports	Documents	Move Equipment	Nx2me	Contact Management	My Account
----------	-----------	----------	---------------	--------------	---------	-----------	----------------	-------	--------------------	------------



Click on a report name to view

Center Inventory: Equipment currently deployed in the field for your center

Unemployed Inventory: Equipment not assigned to a patient

Supplies: Supplies shipped in the last 90 days

QA Release Certification: Equipment approved for use by NxStage Quality Assurance

Patient Information: Names, phone numbers, addresses and status

Patient Therapy: Therapy package, treatment frequency, fluid volume, cartridge type, and suggested dose

Bag Utilization: Tracks supply orders affecting your bag allotment for current and previous months

Nx2me Device Management: Manage patient iPad settings, credentials, and Apple ID requests

Rx Compare: Compare patient prescriptions side by side

Expired Prescriptions: Prescriptions that are expired or nearing expiration

Delivery Schedule: Overview of upcoming patient deliveries

Documents

Documents



Documents

Search

Download

Form	Description
COMM0403	2015 NxStage® Corporate Group Training Programs
FORM112	Prescription Form
FORM236	Patient Status Form
FORM238	Patient Travel Form - Domestic
FORM260	Physician Signature Certification
FORM785	Travel Release Form
FORM816	Patient Travel Form - International

Close

Select an item from the list and click **Download**

- FORM112** – Prescription reference copy
- FORM236** – Notify NxStage of patient hospitalizations, discontinuations, or credit requests
- FORM238** – For U.S. travel of 4 days or more
- FORM260** – 3-year physician signature
- FORM785** – Travel liability release
- FORM816** – Travel outside U.S. after July 10, 2025



Move Equipment

Export to Excel

Report Lost

Move Equipment

View Patient Assets

Designate as Training



Center Inventory

To initiate a move, please select one or more pieces of equipment then click the appropriate action button.

Search **Export to Excel** Move Equipment Designate As Training Report Lost View Patient Assets

Equipment Location	Part	Serial Number	Ship Date	Status	Part Owner

Back

- Select one or more pieces of equipment
- Click **Export to Excel**



Center Inventory

To initiate a move, please select one or more pieces of equipment then click the appropriate action button.

Search

Export to Excel

Move Equipment

Designate As Training

Report Lost

View Patient Assets

Equipment Location	▼ Part	Serial Number	Ship Date	Status	Part Own

Back

- Select one or more pieces of equipment
- Click **Move Equipment**

Move Equipment ← Designate as Training



Center Inventory

To initiate a move, please select one or more pieces of equipment then click the appropriate action button.

[Search](#) [Export to Excel](#) [Move Equipment](#) [Designate As Training](#) [Report Lost](#) [View Patient Assets](#)

Equipment Location	Part	Serial Number	Ship Date	Status	Part Owner

[Back](#)

- Select one or more pieces of equipment
- Click **Designate as Training**



Center Inventory

To initiate a move, please select one or more pieces of equipment then click the appropriate action button.

Search

Export to Excel

Move Equipment

Designate As Training

Report Lost

View Patient Assets

Equipment Location	Part	Serial Number	Ship Date	Status	Part Owner

Back

- Select one or more pieces of equipment
- Click **Report Lost**



Center Inventory

To initiate a move, please select one or more pieces of equipment then click the appropriate action button.

Search

Export to Excel

Move Equipment

Designate As Training

Report Lost

View Patient Assets

Equipment Location

▼

Part

Serial Number

Ship Date

Status

Part Own

Back

- Click **View Patient Assets**
- Review or assign equipment to a patient with no assets assigned

Patient Assets

Please review and assign equipment to new or existing patients where equipment is missing.

Search

Assign Equipment

E-mail Patient

Export to Excel

Export to CSV

1 to 20 of 31

Missing Assets	Patient	Primary Phone	Therapy Type	Cycler SN	Warmer SN	Chassis SN	Control Unit SN	Rx Start Date	Transition To Home Date	Rx Aging
Yes	Frances DemoPatient	(978) 55...	Express (...)					3/31/2015	3/31/2015	2078
Yes	Jacob A DemoPatient	(978) 55...	Express (...)					4/3/2011		3536
Yes	Mark DemoPatient	(978) 55...	Express (...)					9/30/2015		1895



Nx2me

Get Ready to Start

Training

Patients

Order Status

iPad Management



Nx2meLaunchpad
Connected Health

Get Ready To Start | Patients | iPad Management | Training | Order Status

Type of Tr

A

Review Registration

B

Implementation Checklist

Implementation Schedule

Module	Target Complete Date
Module 1: Preparation	10/20/2019
Enroll Patients	10/10/2019
Purchase iPads	10/15/2019

- A. View Center Information
- B. Download Customer Process Checklist

Scroll to view **Implementation Schedule**



 Nx2me Launchpad

Get Ready To Start

Patients

iPad Management

Training

Order Status

Current Nx2me Patients

Patients Awaiting Nx2me Enrollment

Current Nx2me Patients

Anticipated Initial Patients: 2

Search

A Enroll New Patient

B Request ConNx Box Swap

C Nx2me Patient Forms

Export to Excel

Name	DOB	Status	Nx2me Trained?	Enrolled In Nx2me	Package Ship Date
Amy DemoPatient	4/1/1950	Active	Yes	Yes	

- A. Click to select an available patient for enrollment or register a new patient
- B. Select a patient then click to submit a ConNX box swap request for them
- C. Select a patient then click to download their patient form



Nx2me[®]Launchpad
Connected Health

Get Ready To Start

Patients

iPad Management

Training

Order Status

Current Nx2me Patients

Patients Awaiting Nx2me Enrollment

Patients Awaiting Nx2me Enrollment

Anticipated Initial Patient Starts: 2

Search

New Patient

Enroll Selected Patient

Send Connectivity Assessment

Export to Excel

Name	DOB	Status	Email Address	Assessment Sent	Responses Received
Adam HPTaylor	2/5/1975	Active	No Email On Record		
Christopher DemoPatient		Active	demoemail@demo.co	1/11/2024	

- A. Click to register a new patient
- B. Select a patient then click to modify their prescription for Nx2me enrollment
- C. Select a patient then click to send a connectivity assessment

- A.** Copy credentials from an existing iPad
- B.** Guidelines for patient-owned iPads
- C.** Redeploy existing iPad to a new patient



1

1) Preparation

2

Documentation

3

Nurse Implementation Checklist

2

2) iPad Management

3

3) Clinician Portal

4

4) Nx2me App

5

5) Connectivity

6

6) Patient Training & Confirming Home Setup

Docu

ment

Module 1: Preparation

Documentation

Videos

Search

Download

Form

Description

Nurse Implementation Checklist

Checklist of tasks for

Nx2me App Settings Form

This document will

Nx2me Connected Health® Training Materials Rev F

Nx2me Connected

- 1. Select a training **subject**
- 2. Select a training **format**
- 3. Select a training **document** or **video**
- 4. Click **Download** or **Play Video**



Launchpad

Get Ready To Start

Patients

iPad Management

Training

Order Status

Training Order Information

Target Ship Date: 10/5/19

Status: Not Shipped

Clinician Portal Admins

Status: Registered

Portal Access By: 10/5/19

Name: Joe TEST

Name: Test Nurse

Username: joe.test@example.com

Username: nurse@stelmo.com

View your **order status** and Nx2me Supplies that have shipped to your facility within the past 3 months

Contact Management


Creating a Contact

Managing Notifications

Editing a Contact

Managing Permissions

Deactivating a Contact



- Under the **Active Contacts** tab, click **New**
- Enter contact information
- Choose the **Contact type** from the dropdown
- Assign **NxRx Access** if needed
- Click **Submit**

Create Contact

First name

Last name

Email

Phone

Contact type: None

Give NxRx Access: ☐

Submit Cancel

Contact Management ← Editing a Contact



Contact Management

Active Contacts Manage User Permissions Inactive Contacts

Search Manage Notifications New Edit

- Under the **Active Contacts** tab, click **Edit**
- Modify contact information
- Click **Update**

First name Zaid

Last name ADeogun

Email nurse@center.com

Phone (708) 645-1000

Contact type Home Training Nurse

Notification None

Give NxRx access ☐

Update Cancel

Contact Management ← Deactivating a Contact



Contact Management

Active Contacts Manage User Permissions Inactive Contacts

Search Manage Notifications New Edit

- Under the **Active Contacts** tab, click **Edit**
- Change the **Contact type** to:
Inactive DO NOT USE this contact
- Select **Inactive reason**
- Click **Update**

Edit Contact

First name Zaid

Last name ADeogun

Email nurse@center.com


Phone (708) 645-1000

Contact type Inactive DO NOT USE this contact

Inactive reason ☒ Left Company ☐ Unknown ☐ Other

Give NxRx access ☐

Update Cancel




- Under the **Active Contacts** tab, click **Manage Notifications**
- Choose The Notification type from the dropdown
- Click **Update**


Manage Notifications [X]

First name	Olivia
Last name	Gartside
Email	nurse3@center.com
Phone	(517) 787-1893
Notification	None [v]

[✓] Update [✗] Cancel



- Under the **Manage User Permissions** tab, click **Manage Permissions**
- Check the boxes for the desired permissions
- Hover over  to view the permission definition
- Click **Update**





My Account

[Change Password](#)



Get Help	Dashboard	Patients	Welcome Calls	Transactions	Reports	Documents	Move Equipment	Nx2me	Contact Management	My Account
----------	-----------	----------	---------------	--------------	---------	-----------	----------------	-------	--------------------	------------



My Account

User Name

nurse2@center.com

Change Password

Cancel

Change Password

New Password

Confirm Password

Change

Cancel

* Password must contain at least 8 Characters and contain at least: 1 uppercase letter, 1 lowercase letter, 1 number, and 1 special character (!? @\$%)

- Click on **Change Password**
- Enter new password twice to confirm
- Click on **Change**

Thank you

© 2025 Fresenius Medical Care. All Rights Reserved. Fresenius Medical Care, the triangle logo, NxRx and Nx2me are trademarks for which Fresenius Medical Care AG holds rights. All other trademarks are the property of their respective owners. TM1724 Rev A 2025-09